# IVD

Instructions for use (English)

### 1 Intended purpose

The ampliCube Coronavirus Panel is a qualitative in-vitro test to specifically detect the RNA of Middle East Respiratory Syndrome Coronavirus (MERS-CoV), the RNA of novel SARS-CoV-2, SARS-CoV and Bat SARS-like CoV, and the RNA of human coronaviruses (HCoV) (NL63, OC43, 229E, HKU1) in human sputum, swabs, BAL (bronchoalveolar lavage), tracheal secretions or stool.

## 2 Field of application

Middle East Respiratory Syndrome Coronavirus (MERS-CoV) and human coronaviruses (HCoV) including the NL63, OC43, 229E and HKU1 species are part of the Coronaviridae family and are important as the pathogens that cause mild respiratory infections to severe acute respiratory syndrome. About one-third of all colds are caused by HCoV with diarrhoea being a common attendant symptom. MERS-CoV is responsible for Middle East respiratory syndrome, a disease that often manifests as an acute, often severe, respiratory flu-like disease with pneumonia and respiratory distress. Differential diagnostics should consider MERS-CoV particularly if there has been contact with persons infected with MERS-CoV in one of the Middle Eastern countries in the 2 weeks prior to the start of the illness.

SARS coronaviruses such as SARS-CoV-2 (causative pathogen of the COVID-19 pandemic) spread primarily by transmission from person to person via droplets in exhaled air. Symptoms range from fever, cough and dyspnoea to pneumonia and acute respiratory distress syndrome and ultimately death in persons with comorbidities.

The ampliCube Coronavirus Panel is used to detect and differentiate MERS-CoV, SARS-CoV (particularly SARS-CoV-2) and HCoV.

## 3 Test principle

The test is a real-time RT (reverse transcriptase) PCR system. It uses specific primers and labelled probes for transcription of RNA into cDNA and amplification and detection of the RNA of MERS-CoV, SARS-CoV-2, SARS-CoV, Bat SARS-like CoV and HCoV (NL63, OC43, 229E, HKU1).

To ensure that the nucleic acids isolated from the patient sample do not contain any substances that inhibit RT-PCR, an internal control (IC) is added to the sample during the nucleic acid isolation. This IC is transcribed into cDNA, amplified and detected in the same PCR batch. This eliminates false negative test results due to inhibition of the RT-PCR reaction. At the same time, the IC serves as evidence of nucleic acid extraction from the patient sample.

Probes to specifically detect the pathogen-specific nucleic acids are marked with the reporter dyes FAM (MERS-CoV), HEX (SARS-CoV-2, SARS-CoV, Bat SARS-like CoV) and ATTO Rho12 (HCoV) and the probes for detecting the internal control are marked with ATTO 647N. This allows simultaneous detection of all target sequences in a single reaction mixture.

The Ct value (cycle threshold) describes the part of the curve in which the fluorescence rises exponentially above the background value for the first time.

### 4 Reagents

#### 4.1 Package contents

The reagents in one pack are sufficient for 50 assays. Each set of reagents contains:

P&P MIX	<b>150 µI</b> Primer & Probe mix for MERS-CoV, SARS-CoV, HCoV and internal control (green lid)
ENZYME	<b>600 µl</b> enzyme mix ( <b>white lid</b> ) Contains reverse transcriptase and DNA polymerase. ( <b>Component is stained blue</b> .)
CONTROL INT	250 µl internal control (transparent lid)
CONTROL +	170 μl positive control (red lid)
CONTROL -	2 x 1800 µl negative control (blue lid)
INSTRU	1 instructions for use

#### Additionally required reagents, materials and equipment 4.2

- Depending on the used real-time PCR-Cycler, MIKROGEN provides reagents for the colour compensation: MIKROGEN ampliCube Color Compensation (LightCycler® 480 II (Roche) article no. 50502), MIKROGEN ampliCube Color Compensation (cobas z 480 Analyzer (Roche), article no. 50503) or dye calibration kit for CFX96 (Bio-Rad, article no. 50505). If processing at Mic (bms) PCR cycler, MIKROGEN provides Mic assaytemplates.
- Nucleic acid extraction: The following nucleic acid extraction systems are recommended: MagNA Pure® System, Total Nucleic Acid Isolation Kit (Roche) or alphaClean Mag RNA/DNA Kit (MIKRO-GEN) with processing on the M32, M48 or M96 extractor (Biocomma)
- Real-time Cycler: LightCycler® 480 II (Roche), cobas z 480 Analyzer (Roche), CFX96 (Bio-Rad), QuantStudio 5 (Applied Biosystems), Mic (bms), Rotor-Gene Q (Qiagen)
- 96-well PCR plates and films or reaction tubes (PCR-clean): follow the recommendations of the manufacturer of the real-time PCR cycler
- Micropipettes with single-use tips with filter 10 µl, 20 µl, 100 µl and 1000 µl
- Vortex mixer with high rotational speed (recommended 3200 rpm)
- Mini centrifuge
- If necessary, plate centrifuge PBS or H<sub>2</sub>O (*PCR grade*) if using flocked nylon fibre swabs with no transport medium
- Single-use gloves, powder-free •
- Cooling block

#### Shelf life and handling 5

- Store reagents between -25 °C and -18 °C before and after use.
- Repeated thawing and freezing of the components (more than ten times) must be avoided. It is recommended to aliquot the test components after the first thawing.
- Always appropriately chill reagents during the working steps (+2 °C to +8 °C).
- d Keep the kit components away from direct sunlight throughout the test procedure.
- Before starting the test, all reagents must be completely thawed, mixed (briefly vortex) and centrifuged.
- The packages have an expiry date, after which no further guarantee of quality can be given.
- The test must only be performed by trained, authorised and qualified personnel.
- Substantial changes made by the user to the product or the directions for use may compromise the intended purpose of the test specified by MIKROGEN.
- d Cross-contamination can lead to false test results. Add patient samples and controls carefully. Ensure that the reaction mixtures are not transferred to other wells.

#### 6 Warnings and safety precautions

- Only use for in-vitro diagnostics.
- All patient samples must be treated as potentially infectious.
- Suitable single-use gloves must be worn throughout the entire test ø procedure.
- All reagents and materials that come into contact with potentially infectious samples must be treated with suitable disinfectants or be disposed of according to laboratory guidelines. The concentrations and incubation times specified by the manufacturer must be followed.
- Do not replace or mix the reagents with reagents from other kit batches, other MIKROGEN PCR kits or with reagents from other manufacturers.
- Read through and carefully follow all instructions before performing the test. Deviations from the test protocol described in the instructions for use can lead to false results.

## 7 Sampling and preparation of reagents

#### 7.1 Sample material

The starting material for the ampliCube Coronavirus Panel is RNA extracted from sputum, swab, BAL, tracheal secretion or stool of human origin. The quality of the nucleic acid preparation affects the test result. It must be ensured that the extraction method selected is compatible with the real-time PCR technology.

### 7.1.1 Sample preparation

When using a flocked nylon fibre swab without transport medium: 1. Incubate the swab in 0.5 ml PBS (for extraction using the MagNA

- Pure<sup>®</sup> system) or 0.5 ml H<sub>2</sub>O (for extraction using *alpha*Clean Mag RNA/DNA) for 5 min at room temperature.
- 2. Discard the swab and use the PBS or  $H_2O$  suspension for the nucleic acid extraction.

#### 7.2 Extraction of nucleic acids

Extract the nucleic acids from the patient sample and the negative control (NC). We recommend a starting volume for the extraction of 200  $\mu$ l and an elution volume of 50  $\mu$ l or 100  $\mu$ l depending on the extraction system. Extractions (MagNA Pure<sup>®</sup> System (Roche) from 400  $\mu$ l starting material eluted into 100  $\mu$ l showed comparable results. Follow the instructions from the manufacturer of the extraction kit.

1. Thaw the internal control (IC) (transparent lid) and the negative control (NC) (blue lid).

Ensure that the IC and the NC are completely thawed. Mix the IC and the NC before use by briefly vortexing and then briefly centrifuge.

- For the extraction, add 5 µl IC to each patient sample and the NC. The IC should be added to the sample/lysis buffer mix and not directly to the sample material
- 3. Extract the patient sample and the NC. (Note: the NC cannot be used in the PCR without extraction.)
- 4. The positive control is not extracted.

The following nucleic acid extraction system is recommended and was used for the performance assessment:

Extraction system	Sample volume	Elution volume
MagNA Pure <sup>®</sup> 24 (Roche) Total NA Isolation Kit	200 µl	50 µl
M96 Nucleic Acid Extraction System (biocomma) <i>alpha</i> Clean Mag RNA/DNA Kit (MIKROGEN)	200 µl	100 µl

If you would prefer to use other extraction methods, please contact the manufacturer to clarify the compatibility.

#### 7.3 Preparing the master mix

 Thaw the Primer & Probe mix (green lid) and the enzyme mix (white lid). Protect the reagents from light while doing so.
 Ensure that the reagents are completely thawed. Mix the reagents before use by vortexing and then briefly centrifuge.

2. Prepare the master mix using the following pipetting scheme:

Components	Master mix for 1 reaction
Primer & Probe mix	3 µl
Enzyme mix	12 µl
Total volume	15 µl

3. Mix the entire master mix by vortexing and then briefly centrifuge.

4. Produce 15 µl master mix for each PCR reaction.

### 7.4 Preparing the RT-PCR reaction

 Thaw the positive control (PC) (red lid). Ensure that the reagents are completely thawed. Mix the reagents before use by vortexing and then briefly centrifuge.

Components	1 reaction
Master mix from 7.3	15 µl
Sample eluate or eluate of the NC or the PC	10 µl

2. Pipette 10  $\mu$ l of the sample eluate into 15  $\mu$ l master mix.

3. Pipette 10  $\mu I$  of the positive control (not prepared) into 15  $\mu I$  master mix.

4. Pipette 10  $\mu I$  of the eluate of the negative control into 15  $\mu I$  master mix.

Every run must include a positive and a negative control.

Seal the PCR plate with an optically clear adhesive film or seal the reaction tube with the lid provided.



The PCR plates or test tubes must be vortexed for <u>at least 5 sec</u> at maximum speed and then briefly centri-fuged.

PCR test tubes for the <u>Mic PCR cycler</u> must be vortexed for at least 10 sec at maximum speed.

### 8 Programming the real-time cycler

The *ampli*Cube Coronavirus Panel was evaluated with the LightCycler<sup>®</sup> 480 Instrument II (Roche) and validated at cobas z 480 Analyzer (Roche), CFX96 (Bio-Rad), QuantStudio 5 (Applied Biosystems), Mic (bms) and Rotor-Gene Q (Qiagen).

### 8.1 Setting the detection channels

#### LightCycler® 480 Instrument II (Roche)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red
Reporter dye	FAM	HEX	ATTO Rho12	ATTO 647N
Excitation	465 nm	533 nm	533 nm	618 nm
Emission	510 nm	580 nm	610 nm	660 nm

For the LightCycler<sup>®</sup> 480 II, it is necessary to first use the Color Compensation kit (article no. 50502), which is provided by MIKROGEN.

#### cobas z 480 Analyzer (Roche)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red
Reporter dye	FAM	HEX	ATTO Rho12	ATTO 647N
Excitation	465 nm	540 nm	540 nm	610 nm
Emission	510 nm	580 nm	610 nm	670 nm

For the cobas z 480 Analyzer (Roche), it is necessary to first use Color Compensation (article no. 50503), which is provided by MIKROGEN.

#### CFX96 (Bio-Rad)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red
Reporter dye	FAM	HEX	ATTO Rho12	ATTO 647N
Mode	all channels			

The data acquisition and evaluation at CFX96 is carried out using the method *calc / data acquisition mode: all channels*. A dye calibration of the CFX96 (Bio-Rad) must be carried out in advance for ATTO Rho12 and ATTO 647N, which is provided by MIKROGEN (CFX96 (Bio-Rad) dye calibration set, article no. 50505).

#### QuantStudio 5 (Applied Biosystems)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red
Reporter dye	FAM	HEX	ATTO Rho12	ATTO 647N
Excitation	X1 / 470 nm	X2 / 520 nm	X4 / 580 nm	X5 / 640 nm
Emission	M1 / 520 nm	M2 / 558 nm	M4 / 623 nm	M5 / 682 nm
Quencher	[none]	[none]	[none]	[none]

Choose under settings 1. Run mode "standard", 2. Reference dye "none", 3. Experiment type "custom".

#### Mic (bms)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red

MIKROGEN provides you with validated Mic assay-templates for processing. Please use exclusively MIKROGEN Mic assay-templates. The current MIKROGEN Mic-Assay-Templates are freely available for download on the MIKROGEN homepage (www.mikrogen.de).

#### Rotor-Gene Q (Qiagen)

	MERS-CoV	SARS-CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colour	Green	Yellow	Orange	Red
Reporter dye	FAM	HEX	ATTO Rho12	ATTO 647N
Excitation	470 nm	530 nm	585 nm	625 nm
Emission	510 nm	555 nm	610 nm	660 nm
Gain	5	5	5	5

### 8.2 PCR program

Reverse transcription	50 °C	8 minutes
Denaturation	95 °C	3 minutes
Amplification	45 cycles	
Denaturation	95 °C	10 seconds
Annealing/elongation	60 °C	45 seconds

Essential information about programming the different real-time cycles can be found in the instructions for the cycler used. For specific information about programming the real-time PCR cycler when using the *ampli*Cube Coronavirus Panel, please contact the manufacturer.

### 9 Results

The data analysis on the LightCycler<sup>®</sup> 480 II is carried out using the *Abs Quant/2nd Derivative Max* method.

#### 9.1 Validation

- 1. The negative control must be below the threshold. The internal control (IC) must have a positive curve for the negative control. If the negative control has a positive curve (contamination) or if the IC is not valid in the negative control, the test run cannot be analysed.
- The positive control must have a positive curve. The Ct value for the positive control must be < 33. A positive control outside this range indicates that there is a problem with the amplification.
- 3. The internal control (IC) must have a positive curve for negative samples. The signal for the IC for a patient sample must be compared to the signal of the IC in the extracted negative control. A difference of >+3 for the Ct value of the IC of a sample compared to the IC of the negative control or the absence of an IC signal for the sample may indicate significant inhibition of the RT-PCR reaction. In these cases a negative test result is invalid.

#### 9.2 Evaluation

The data can be analysed with the corresponding PCR cycler software or a software solution for automated PCR analysis and interpretation specifically supported by MIKROGEN. Using a LightCycler<sup>®</sup> 480 II the analysis can be carried out either using the *Abs Quant/2nd Derivative Max* method (recommended) or the *Abs Quant/Fit Points* method. Additional information and corresponding instructions are available from MIKROGEN upon request.

		MERS-CoV	SARS- CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Colo	our	Green	Yellow	Orange	Red
Rep	orter dye	FAM	HEX	ATTO Rho12	ATTO 647N
h	LightCycler 480 II	465-510	533-580	533-610	618-660
with	cobas z 480	465-510	540-580	540-610	610-670
uo	CFX96	FAM	HEX	ROX	Cy5
uati	QuantStudio 5	FAM	VIC	ROX	ATTO 647N
Evaluation	Mic	FAM	HEX	ROX	Cy5
Ē	Rotor-Gene Q	green	yellow	orange	red

Signals above the threshold are evaluated as positive results. Empty fields in the table are considered negative results.

Colour	MERS-CoV	SARS- CoV-2 SARS-CoV Bat SARS- like CoV	HCoV	Internal control (IC)
Green	Positive			
Yellow		Positive		
Orange			Positive	
Red				Positive*

\* In the case of positive signals in the detection channels for the pathogens, the signal for the internal control is not required for the test interpretation. A high pathogen load in the patient sample can lead to a reduced or missing signal for the internal control.

### 10 Limits of the method, restrictions

- Test results must always be viewed in the context of the clinical situation. Therapeutic consequences of the findings must be linked to the clinical data.
- A negative MERS-CoV, SARS-CoV-2/SARS-CoV / Bat SARS-like CoV or HCoV test result cannot rule out infection with the particular pathogens.

### **11 Performance characteristics**

11.1 Diagnostic sensitivity and specificity

The sensitivity and specificity were determined using defined positive and defined negative samples.

#### Table 1: Defined positive samples

ampliCube	MERS-CoV	SARS-CoV-2	HCoV
Coronavirus Panel			
	(n=4)	(n=126)	(n=17)
Negative	0	3	0
Positive	4	123	17
Sensitivity [%]	100	97.62	100
<b>CI</b> [%] Confidence interval	51.01 – 100	93.23 - 99.19	81.57 – 100

Table 2: Defined negative samples

ampliCube	MERS-CoV	SARS-CoV-2	HCoV
Coronavirus Panel			
	(n=24)	(n=490)	(n=24)
Negative	24	486	24
Positive	0	4	0
Specificity [%]	100	99.18	100
CI [%] Confidence interval	86.20 – 100	97.92 - 99.68	86.20 – 100

#### 11.2 Analytical sensitivity

11.2.1 Limit of Detection (LoD) in copies/PCR

The limit of detection (LoD) of the *ampli*Cube Coronavirus Panel was determined using a dilution series of gBlock DNA of known concentration in a LightCycler<sup>®</sup> 480 II System (Roche). The 95% limit of detection was determined using probit regression analysis with CombiStats<sup>™</sup> Version 5.0 software (Council of Europe). The limit of detection was then confirmed with 20 replicates at the LoD.



### *ampli*Cube Coronavirus Panel Instructions for use (English)

Table 3a: Limit of detection (LoD) for MERS and HCoV in copies/PCR

	MERS-CoV	HCoV
LoD [copies/PCR] 95% limit of detection	4.14	12.41*
95% CI [copies/PCR] confidence interval	2.75 – 9.43	8.93 – 22.59

\* Information refers to HCoV OC43; HCoV 229E, NL63 and HKU1 were also tested. The limits of detection are between 20.54 and 89.48.

Table 2b. Limit of datastion (		for SABS Cold		) in conice/DCD	
Table 3b: Limit of detection (	LOD	101 SARS-COV 1	LOD	) in copies/PCR	

	SARS-CoV-2	SARS-CoV	Bat SARS-like CoV
LoD [copies/PCR] 95% limit of detection	2.9	6.74	6.35
95% CI [copies/PCR] confidence interval	1.9 - 6.4	4.12 – 14.68	4.14 - 14.43

11.2.2 Limit of Detection (LoD) SARS-CoV-2 in International Units In addition, the LoD for the SARS-CoV-2 target regions was determined with three independent dilution series (7 dilution levels each) using the WHO standard (First WHO International Standard for SARS-CoV-2 RNA (NIBSC code: 20/146)) including nucleic acid extraction (*alpha*Clean Mag RNA/DNA Kit on M96 Biocomma) on a LightCycler<sup>®</sup> 480 II system (Roche) and subsequent probit regression analysis (CombiStats<sup>™</sup> version 6.0 software (Council of Europe)). The detection limit of the *ampl*/Cube Coronavirus Panel is given in International Units (IU) per µl WHO standard.

Tabelle 3c: Limit of detection (LoD) for SARS-CoV-2 in IU/µI WHO standard

	SARS-CoV-2
LoD [IU/µI] 95% limit of detection	2.72
<b>95% CI</b> [IU/µI] confidence interval	1.53 – 6.55

#### 11.3 Analytical specificity

The BLAST trial (www.ncbi.nlm.nih.gov/blast/) shows that the selected primers and probes of the *ampli*Cube Coronavirus Panel specifically detect the selected pathogens.

Furthermore, the specificity was determined by studying genomic DNA/RNA of other human pathogenic bacteria and viruses.

Table 4: Bacteria and viruses that were tested in order to demonstrate the analytical specificity of the *ampli*Cube Coronavirus Panel.

Bacteria	Viruses	
Bordetella holmesii	Adenovirus A	
Bordetella parapertussis	Cytomegalovirus	
Bordetella pertussis	Enterovirus	
Chlamydia pneumoniae	Epstein-Barr Virus	
Mycoplasma pneumoniae	Influenza A	
	Influenza B	
	Measles virus	
	Mumps virus	

Measles virus
Mumps virus
Parainfluenzavirus
Respiratory syncytial virus (RSV) A
Rhinovirus 58 (A)

None of these samples generated a positive signal. The primers and probes used in the *ampli*Cube Coronavirus Panel showed no cross-reactions with the pathogens listed in Table 4. The internal control (IC) was valid in all tests.

#### 11.4 Equivalence of different sample material

The coefficient of variation (CV) was determined for the Ct value between water and the extract of the particular sample material after the addition of plasmid DNA of known concentration.

 Table 5: Equivalence of different sample material

	MERS- CoV	SARS- CoV-2	SARS- CoV	Bat SARS- like CoV	HCoV
CV [%] (BAL, H <sub>2</sub> O)	1.37	2.11	1.62	2.21	0.58
CV [%] (sputum, H <sub>2</sub> O)	2.42	2.32	1.54	3.42	1.91
CV [%] (swab, H <sub>2</sub> O)	1.51	1.28	1.47	1.50	1.08
CV [%] (stool, H <sub>2</sub> O)	1.94	1.33	1.28	1.98	0.85



The coefficient of variation (CV), based on the Ct value (cycle threshold) between water and the nucleic acid extracts (obtained for the various sample materials), was  $\leq 3.42\%$  for all target genes.

#### **12 Literature**

- L.-Y. Chang et. al. (2006): Lack of Association between Infection with a Novel Human Coronavirus (HCoV), HCoV-NH, and Kawasaki Disease in Taiwan. The Journal of Infectious Diseases, Jan. 2006, 193:283-6, 0022-1899/2006/19302-0015\$15.00
- L. Louie et. al. (2006): Detection of Severe Acute Respiratory Syndrome Coronavirus in Stool Specimens by Commercially Available Real-Time Reverse Transcriptase PCR Assays. Journal of Clinical Microbiology, Nov. 2006, Vol. 44, No. 11, p. 4193-4193, doi:10.1128/JCM.01202-06
- F. Esper et. al. (2005): Evidence of a Novel Human Coronavirus That Is Associated with Respiratory Tract Disease in Infants and Young Children. The Journal of Infectious Diseases, Feb. 2005, 191:492-8, 0022-1899/2005/19104-0002\$15.00
- F. Esper et. al. (2010): Human coronaviruses are uncommon in patients with gastrointestinal illness. Journal of Clinical Virology, March 2010, Vol. 48, p. 131-133, doi:10.1016/j.jcv.2010.03.007
- M. Jevšnik et. al. (2013): Detection of human coronaviruses in simultaneously collected stool samples and nasopharyngeal swabs from hospitalized children with acute gastroenteritis. Virology Journal, 2013, 10:46, doi:10.1186/1743-422X-10-46
- I. M. Mackay et. al. (2015): MERS coronavirus: diagnostics, epidemiology and transmission. Virology Journal, Dec. 2015, 12:222, DOI 10.1186/s12985-015-0439-5
- M. C. Nunes et. al. (2014): Clinical Epidemiology of Bocavirus, Rhinovirus, Two Polyomaviruses and Four Coronaviruses in HIV-Infected and HIV-Uninfected South African Children. Plos One, Feb. 2014, Vol. 9, Issue 2, e86448
- L. A. Sipulwa (2016): Molecular characterization of human coronaviruses and their circulation dynamics in Kenya, 2009–2012. Virology Journal, 13:18, DOI 10.1186/s12985-016-0474-x
- 9. I. Wilhelmi (2003): Viruses causing gastroenteritis. Clinical Microbiology and Infection, April 2003, Vol. 9, No. 4, p. 247-262
- 10. F. Wu et. al. (2020): A new coronavirus associated with human respiratory disease in China. Nature (2020) DOI:10.1038/s41586-020-2008-3
- R. Ralph et. al. (2020): 2019-nCoV (Wuhan virus), a novel Coronavirus: human-to-human transmission, travel-related cases, and vaccine readiness. J Infect Dev Ctries 2020; 14(1):3-17
- W. G. Carlos et. al. (2020): Novel Wuhan (2019-nCoV) Coronavirus. Am J Respir Crit Care Med 2020 Vol. 201, P7-P8
- S.-Q. Deng et H.-J. Peng (2020): Characteristics of and Public Health Responses to the Coronavirus Disease 2019 Outbreak in China. J. Clin. Med. 2020, 9, 575-584
- Muenchhoff Maximilian, Mairhofer Helga, Nitschko Hans, Grzimek-Koschewa Natascha, Hoffmann Dieter, Berger Annemarie, Rabenau Holger, Widera Marek, Ackermann Nikolaus, Konrad Regina, Zange Sabine, Graf Alexander, Krebs Stefan, Blum Helmut, Sing Andreas, Liebl Bernhard, Wölfel Roman, Ciesek Sandra, Drosten Christian, Protzer Ulrike, Boehm Stephan, Keppler Oliver T. Multicentre comparison of quantitative PCR-based assays to detect SARS-CoV-2, Germany, March 2020. Euro Surveill. 2020;25(24):pii=2001057.https://doi.org/10.2807/15607917.ES.2020.25.24.2 001057
- 15. Ute Eberle, Clara Wimmer, Ingrid Huber, Antonie Neubauer-Juric, Giuseppe Valenza, Nikolaus Ackermann, Andreas Sing (for the Bavarian SARS-CoV-2-Public Health Laboratory Team), Comparison of nine different commercially available molecular assays for detection of SARS-CoV-2 RNA, European Journal of Clinical Microbiology & Infectious Diseases, Springer, published online: 29 January 2021, https://doi.org/10.1007/s10096-021-04159-9

We will be pleased to send you additional literature on request.

## **13 Explanation of symbols**

$\Sigma$	Content is sufficient for <n> formulations Number of formulations</n>				
P&P MIX	Primer & Probe mix				
ENZYME	Enzyme mix				
CONTROL INT	Internal control				
CONTROL +	Positive control				
CONTROL -	Negative control				
INSTRU	Instructions for use				
	Follow the instructions for use				
CONT	Contents, contains				
IVD	In-vitro diagnostic medical device				
LOT	Batch number				
REF	Order number				
2	Use by Expiry date				
x°C y°C	Store between x°C and y°C				
	Manufacturer				

## 14 Manufacturer and version data

ampliCube	e Coronavirus Panel		Article no. 50142
Instruction Valid from	ns for use		GAACCV003EN 2022-05
	MIKROGEN GmbH Floriansbogen 2-4 82061 Neuried Germany Tel. Fax Email Internet	+49 89 54801-0 +49 89 54801-100 mikrogen@mikrog www.mikrogen.de	gen.de
			CE



