

XCeloSeq® Sarcoma Fusion Kit

SEQ014

Product Description

The XCeloSeq Sarcoma Fusion Kit contains a pool of targeted RNA enrichment primers located in conserved fusion partners for identification of both known and unknown fusions from RNA. These primers are designed for use only with XCeloSeq Targeted RNA Core Reagents (GF031). Together they allow for the generation of high quality, high-complexity next-generation sequencing libraries that are suitable for use with Illumina® next-generation sequencing instruments.

Kit Contents

Component	Tube Colour	Cap Colour	Storage	Part Code
Sarcoma Fusion Kit – Outer Pool	Transparent	Orange	-20°C	PC0250
Sarcoma Fusion Kit – Inner Pool	Transparent	Black	-20°C	PC0251

Specifications

Gene Targets	26
Targeting Primers [%]	146
Recommended Input Quantity [*]	5-200 ng total FFPE-RNA
	5-100 ng high quality RNA
Recommended Reads Per Sample	1,500,000
Hands on Time	2.0 Hours
Total Protocol Time	7.25 hours

[%]An additional 8 QC primers are included

^{*}Higher quantities within this range will improve maximum sensitivity. The product supports capture with down to 1.0 ng of RNA, however this is not recommended as it will lead to reduced sensitivity. Cell-free RNA and total cell-free nucleic acids may be used as alternative starting materials, however fusion detection sensitivity will be lower due to cell-free RNA concentrations typically being very low, when using this material maximising starting input quantity will help ensure the best possible results.

Assay Targets

Gene	Accession	Exon(s)	Fusion Direction
ALK	NM_004304	19, 20, 21, 22	5'
CAMTA1	NM_015215	8, 9, 10	5'
CCNB3	NM_033031	2, 3, 4, 5, 6	5'
CIC	NM_015125	19, 20	3'
EPC1	NM_025209	9, 10, 11	3'
EWSR1	NM_005243	4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14	3'
FOXO1	NM_002015	1, 2, 3	5'
		1, 2, 3	3'
FUS	NM_004960	4, 5, 6, 7, 8, 9, 10, 11, 14	3'
GLI1	GLI1	4, 5, 6, 7	5'
		4, 5, 6, 7	3'
HMGA2	NM_003483	1, 2, 3, 4, 5	3'
JAZF1	NM_175061	2, 3, 4	3'
MEAF6	NM_001270875	4, 5	3'
MRTFB	NM_014048	11, 12, 13	5'
NCOA2	NM_006540	11, 12, 13, 14	5'
NTRK3	NM_002530	4, 7, 10, 12, 13, 14, 15, 16	5'
		13, 14, 15	3'
	NM_001007156	15	5'
PDGFB	NM_002608	2, 3	5'
PLAG1	NM_002655	1, 2, 3, 4	5'
ROS1	NM_002944	31, 32, 33, 34, 35, 36, 37	5'
SS18	NM_001007559	10, 11	5'
		4, 5, 6, 8, 9, 10	3'
STAT6	NM_001178078	1, 2, 3, 4, 5, 6, 7, 16, 17, 18, 19	5'
TAF15	NM_139215	6, 7	5'
		5, 6, 7	3'
TCF12	NM_207036	4, 5, 6	3'
TFE3	NM_006521	3, 4, 5, 6	5'
TFG	NM_006070	4, 5, 6, 7	5'
		6	3'
USP6	NM_004505	1, 2, 3	5'
YWHAE	NM_006761	5	3'

Additional Information

Please refer to “XCelSeq Targeted RNA Enrichment Protocol” for instructions for use.

Limitations of Use

For Research Use Only (RUO)

This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals. SDS sheets relevant to this product are available upon request.

Illumina® is a registered trademark of Illumina, Inc.

Customer Service & Sales Enquiries
Telephone: +44 (0)1865 407 400
Email: sales@genefirst.com
Document: INS1401v1.0

GeneFirst Limited
Building E5, Culham Science Centre,
Abingdon, Oxfordshire. OX14 3DB UK

© 2021 GeneFirst, Ltd. All rights reserved.
www.genefirst.com
Page 3