



XPERT® BLADDER CANCER MONITOR*

Test Reagent Kit	Xpert® Bladder Cancer Monitor
Test Kit Catalog Number	GXBLAD-CM-CE-10
Technology	Real-time PCR
Targets	Five mRNA targets UPK1B, IGF2, CRH, ANXA10, ABL1
Batch or On-Demand	On-Demand
Minimum Batch Size	1
Sample Type	Voided Urine
Sample Extraction	Automated/Integrated
TAT	90 minutes
Hands-On Time	<2 minutes
Controls: Sample Adequacy Control	ABL1
Controls: Probe Function/Detection	Probe Check Control
Controls: Internal Control	Armored RNA

Xpert Bladder Cancer Monitor Performance vs. Cystoscopy/Histology

	(N=257)	Overall	Low Grade	High Grade
Sensitivity		75.0%	63.2%	84.0%
Specificity		80.6%	N/A	N/A

Negative Predictive Value	93.9%
Positive Predictive Value	44.6%
Sample Storage & Stability	Urine samples in Xpert Urine Transport Reagent tubes are stable up to 7 days at 2–28 °C
Kit Storage	2–28° C
Additional Required Materials	Xpert Urine Transport Reagent Kit <i>Catalog Number:</i> GXUTR-CE-30
Commercial Controls	Xpert Bladder Cancer QC Panel from Molecular Quality Controls, Inc (MMQC) <i>Catalog Number:</i> C104 <i>Contains:</i> 6 x 4mLs of Xpert Bladder Cancer low positive control 6 x 4mLs of Xpert Bladder Cancer negative control



Not for distribution in the U.S.

CORPORATE HEADQUARTERS904 Caribbean Drive
Sunnyvale, CA 94089 USATOLL FREE 1.888.336.2743
PHONE 1.408.541.4191
FAX 1.408.541.4192**EUROPEAN HEADQUARTERS**Vira Solelh
81470 Maurens-Scopont FrancePHONE 33.563.82.53.00
FAX 33.563.82.53.01
EMAIL cepheid@cepheideurope.fr

www.Cepheidinternational.com

