1. Intended Use

For *In Vitro* Diagnostic Use.

The IDS-iSYS Insulin-like Growth Factor-I (IGF-I) assay is an *in vitro* diagnostic device intended for the quantitative determination of IGF-I in human serum or plasma on the IDS-iSYS Multi-Discipline Automated System. Results are to be used in conjunction with other clinical and laboratory data to assist the clinician in the assessment of growth disorders.

2. Summary and Explanation

Insulin-like growth factor-I (IGF-I) is a polypeptide of 70 amino acids (7650 Daltons), and is one of a number of related insulin like growth factors present in the circulation. The molecule shows approximately 50% sequence homology with proinsulin and has a number of biological activities similar to insulin. The peptide is growth hormone (GH) dependent to a high degree, but there is growing evidence of GH-independent secretion. IGF-I has numerous growth-promoting effects, including mitogenic effects and the promotion of cartilage sulphation. It also mediates growth promoting actions of growth hormone on skeletal and other body tissues.

Almost all (>95%) of serum IGF-I circulates bound to specific IGF binding proteins, of which six classes (IGFBPs 1-6) are now recognised¹. IGFBP-3 is thought to be the major binding protein of IGF-I, forming a ternary complex of 140 000 Daltons with IGF-I and an acid labile subunit².

The measurement of serum IGF-I is of recognised value in children with growth disorders and in the diagnosis and monitoring of acromegaly³. IGF-I concentrations change with age, nutritional status, body composition and GH secretion.

A single basal IGF-I determination is useful in the assessment of short stature in children and in nutritional support studies of acutely ill patients. For the diagnosis of acromegaly, a single IGF-I determination is considered more reliable than a random GH determination⁴.

3. Method Description

The assay is based on chemiluminescence technology. 10 µL of patient sample are incubated with an acidic solution to dissociate IGF-I from the binding proteins. A portion of this, along with neutralisation buffer is incubated with a biotinylated anti-IGF-I monoclonal antibody, and an acridinium labelled anti-IGF-I monoclonal antibody. Streptavidin labelled magnetic particles are then added. The magnetic particles are captured using a magnet and a wash step performed to remove any unbound analyte. Trigger reagents are added; the resulting light emitted by the acridinium label is directly proportional to the concentration of IGF-I in the original sample.

4. Warnings and Precautions

The IDS-iSYS IGF-I assay is for *in vitro* diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in these Instructions For Use (IFU). Immunodiagnostic Systems Limited (IDS) will not be held responsible for any loss or damage (except as required by statute), howsoever caused, arising out of non-compliance with the instructions provided.

CAUTION: This kit contains material of animal origin. Handle kit reagents as if capable of transmitting an infectious agent. Appropriate precautions and good laboratory practice must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Human materials

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV I and II), Hepatitis B surface antigen, antibody to Hepatitis C, and found negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled according to Biosafety Level 2.

Reagents containing Sodium Azide

Some reagents in this kit contain sodium azide (NaN_3) which may react with lead, copper or brass plumbing to form highly explosive metal azides. On disposal, flush with large volumes of water to prevent azide build up. Calibrators [CAL] contain $\text{NaN}_3 >0.1\%$ (w/w) (<8%).

Classification under CLP:

Acute Tox. 4

Aquatic Chronic 3

Hazard statements:

- EUH032 Contact with acids liberates very toxic gas
- H302 Harmful if swallowed
- H412 Harmful to aquatic life with long lasting effects

Precautionary statements:

- P264 Wash hands thoroughly after handling.
- P270 Do not eat, drink or smoke when using this product.
- P273 Avoid release to the environment.
- P301/310 IF SWALLOWED: Immediately call a POISON CENTER or doctor.
- P330 Rinse mouth.
- P501 Dispose of contents/container to hazardous or

special waste collection point.

5. Handling Precautions

Apart from the calibrators, which are lyophilised, the reagents provided in the kit are ready to use. Refer to the Calibrator section of the procedure for reconstitution methodology.

Before a new cartridge is loaded onto the IDS-iSYS Multi-Discipline Automated System (System), the magnetic particle container requires mixing by the operator with a brisk rotation between the palms of the hands in a back and forth motion or via the IDS-iSYS XPrep device. This will re-suspend the magnetic particles that have settled during shipment. Care must be taken to avoid foam formation in the reagent cartridge.

Shelf Life and Storage of Reagents

Store the cartridge and calibrators in an **upright** position in the dark at 2 - 8 °C. **Do not freeze** the cartridge or the calibrators.

Reagent shelf life	Cartridge	Calibrators
Before opening at 2 - 8 °C	To the expiry date	
Cartridge - After opening at 2 - 8 °C	28 Days	N/A
Calibrators - After reconstitution at -20 °C in the original vials	N/A	7 weeks
On board the System *	10 Days	3 Hours

* Continuous on board stability

6. Sample Collection and Storage

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin and potassium EDTA) samples. Samples should be separated as soon as possible after collection. Store samples at -20 °C or below. Avoid repeated freeze - thaw of samples.

Note:

- Some commercially available sample collection tubes may affect the results of testing in particular cases.
- Follow the blood collection tube manufacturer's recommendations for handling and processing the samples.
- Samples containing particulate matter must be centrifuged before performing the assay. Centrifuged samples with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified samples without the lipemic material
- Before performing assays, make sure that samples, calibrators and controls are at room temperature (18 - 22 °C).
- To minimise possible evaporation effects, samples, calibrators and controls should be measured within 3 hours of being placed on the System.
- Do not use heat-inactivated samples.
- Each laboratory should follow the guidelines or requirements of local, state, and/or federal regulations or accrediting organizations to establish its own specimens handling and storage stability. For guidance on appropriate practices, please refer to the CLSI GP44-A4, "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests".

7. Materials

Materials Provided

Reagent Cartridge

MPT1	Magnetic particles coated with Streptavidin in a phosphate buffer containing sodium azide as preservative (<0.1%), 1 bottle, 2.7 mL
CONJ	Anti-IGF-I labelled with an acridinium ester derivative, in buffer containing bovine proteins with sodium azide as preservative (<0.1%), 1 bottle, 4.25 mL
Ab-BIOT	Anti-IGF-I labelled with biotin, in buffer containing bovine and mouse proteins with sodium azide as preservative (<0.1%), 1 bottle, 4.25 mL
BUF	A phosphate buffer containing Tris and IGF-II, with sodium azide as preservative (<0.1%), 1 bottle, 12.5 mL
DIL	An acidic sodium chloride solution. 1 bottle, 40 mL

Calibrators

CAL A	Bovine serum albumin in a buffer matrix containing IGF-I and sodium azide as preservative (<8%), 1 each of 2 concentration levels, 1.0 mL per bottle
CAL B	

Mini CD

Contains IFU for IDS-iSYS reagents, CRY files and Certificate of Analysis.

Materials Required But Not Provided

Standard Materials

IS-310400	IDS-iSYS Multi-Discipline Automated System
IS-3930	IDS-iSYS Insulin-like Growth Factor-I (IGF-I) Control Set
IS-CC100	IDS-iSYS Cuvettes
IS-CS100	IDS-iSYS System Liquid (Syst. I)
IS-CW100	IDS-iSYS Wash Solution (Wash S)
IS-CT100	IDS-iSYS Triggers Set

Optional equipment / materials

IS-3935	IDS-iSYS iSYS Insulin-like Growth Factor-I IGF-I Calibration Verifiers
IS-XP01	IDS-iSYS XPrep
	Precision Pipetting Devices

IS-6010	IDS-iSYS Cartridge Check System (CCS)
IS-CSC105	Sample Cups (500 µL)
IS-DS100	D-SORB solution

8. Assay Procedure

Reagent Cartridge

The reagents provided in the cartridge are ready to use, with the exception of the calibrators. Before a new cartridge is loaded on board the System, mix the magnetic particles container with a brisk rotation between the palms of the hands in a back and forth motion or via the IDS-iSYS XPrep device. Avoid foam formation.

The barcode is read when the cartridge is loaded on the reagent tray. If the label cannot be read by the System barcode reader, a manual procedure exists to enter the barcode data (see the IDS-iSYS User Manual).

The System automatically performs mixing of the magnetic particles for all reagent cartridges loaded on-board to maintain homogeneity. The cartridge should be loaded into the reagent tray at least 40 minutes prior to starting the assay. If the cartridge is removed from the reagent tray, store the cartridge upright at 2 - 8 °C in the dark.

Do not perform the cartridge level scan on the IDS-iSYS IGF-I cartridge as this will result in less than 100 tests being obtained from the cartridge.

Calibrators

The calibrators of this kit are traceable to the WHO International Standard for IGF-I, code 02/254.

The IGF-I calibrators are lyophilised. Reconstitute immediately before use. Add 1.0 mL of distilled or deionised water to each bottle. Replace the stopper. Leave for 10 minutes to reconstitute with occasional gentle mixing by hand. Care should be taken to avoid the formation of foam.

Calibrators can be loaded onto the System using either the sample cups or calibrator/control racks (30/31/32).

Calibrator/control racks: place the calibrator vials into the appropriate rack and insert onto the System and proceed according to the instructions of the System User Manual. Discard the material in the calibrator vials after use. DO NOT re-use calibrator vials that have been inserted onto the System.

Sample cups: Pipette approximately 200 µL of calibrators into sample cups and place onto the system. Proceed according to the instructions of the System User Manual. Discard the material in the sample cups after use. DO NOT return material to the calibrator vial.

If calibrators are to be used more than once, calibrators should be stored at -20 °C or lower within 15 minutes in their original vials following reconstitution. When re-using the stored calibrator vials, bring the calibrators to room temperature mix well. Ensure that the Calibrators are at room temperature before they are placed on the System. Calibrators should be placed on the System within 30 minutes of thawing. Thawed Calibrators should not be re-frozen.

Calibration Frequency

A new calibration is required:

- Upon loading each new lot of cartridges.
- When the control values do not fall within the defined ranges.
- When the calibration interval of 10 days has expired.
- After System service.

Verification of the calibration is automatic and managed by the System.

System Calibration

The two IGF-I calibrators are required to perform the adjustment of the master curve. The calibrators are supplied with the kit and calibrators from another lot must not be used. All levels of IDS-iSYS IGF-I Control Set (IS-3930) MUST be measured at the same time as the calibrators to validate the master curve adjustment.

All data required for the calibration of the cartridge can be found on the mini CD. Use calibrator levels A and B to adjust the master curve to the reagents on board the System. Check for the presence of an IGF-I cartridge on the reagent tray and the availability of the cartridge master curve in the database. If the data for the lot of calibrators is not available on board the System, load the data using the mini CD provided with the reagent kit.

Start the immunoassay calibration on the System according to the System User Manual. The calibration is carried out using 3 replicates of each level. One replicate may be removed to meet the calibration requirements. All control levels must also be measured in duplicate to calibrate the assay. Verify and approve the calibration according to the calibration status displayed in the calibration window and discard the remaining volume of calibrators and controls from the sample tray after use.

9. Quality Control

Use the IDS-iSYS IGF-I Control Set for quality control. To ensure validity of results at least three controls with varying levels of IGF-I should be measured. Other suitable control material can be used in addition to the IDS-iSYS IGF-I Control Set. Controls should be tested at (or near) the beginning of every run containing patient samples and also during every calibration in accordance with local, state and/or federal regulations or accreditation requirements and your laboratory's quality procedure. Refer to the IDS-iSYS IGF-I Control Set IFU for preparation and handling instructions.

Determination of Sample IGF-I levels

Process samples according to the System User Manual.

10. Calculation of Results

The concentration of IGF-I within each sample is calculated automatically. The display of the concentrations (screen or printed) is produced upon user request.

The IDS-iSYS IGF-I assay uses a 4-parameter logistic curve fit (4PL) to calculate the IGF-I concentrations.

11. Measurement Range (Reportable Range)

The reportable range of the assay is 10 to 1200 ng/mL. Any value that reads below 10 ng/mL should be reported as "*< 10 ng/mL*".

12. Dilution

Samples with IGF-I concentrations above the reportable range should be diluted manually with low concentration human serum sample in a ratio of 1 in 2. The results for diluted samples must be multiplied by the dilution factor 2 and corrected for the concentration of the low sample.

13. Limitations of Use

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The hook effect was tested using concentrations of IGF-I up to 20000 ng/mL. No hook effect was observed.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁶. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The following substances do not interfere in the IDS-iSYS IGF-I assay when the concentrations presented in the following table are below the stated threshold.

Potentially Interfering Agent	Threshold Concentration
Lipid	3000 mg/dL
Bilirubin	20 mg/dL
Haemoglobin	500 mg/dL
Biotin	300 nmol/L
IGFBP1	5000 ng/mL
IGFBP2	5000 ng/mL
IGFBP3	20000 ng/mL
IGFBP4	5000 ng/mL
IGFBP5	5000 ng/mL
IGFBP6	5000 ng/mL

14. Expected Values

The following ranges were determined using the IDS-iSYS IGF-I assay in a large, multicentre (n=12) cohort (n=15014) from the United States, Canada and Europe⁶. Even though in order to best reflect local conditions it is recommended that each laboratory establish its own expected reference range based upon its own patient population, the following ranges should be considered as appropriate normative data for this assay. Laboratories who wish to have more detailed ranges with data in smaller increments are referred to the online Supplemental data in the relevant publication⁶.

IGF-I Males						IGF-I Females					
LMS Parameters ($\lambda = 0.4568$)			Percentiles			LMS Parameters ($\lambda = 0.4595$)			Percentiles		
Age, y	μ	σ	2.5%	50%	97.5%	Age, y	μ	σ	2.5%	50%	97.5%
0	77.3625	0.4262	27.0	77.4	157.0	0	58.5742	0.4662	17.9	58.6	125.6
1	83.0813	0.4187	29.7	83.1	166.8	1	62.3128	0.4589	19.5	62.3	132.3
2	92.5995	0.4113	33.9	92.6	183.9	2	69.1965	0.4516	22.2	69.2	145.4
3	104.0358	0.4038	39.0	104.0	204.5	3	78.9225	0.4444	25.9	78.9	164.2
4	115.7263	0.3964	44.3	115.7	225.0	4	91.1878	0.4371	30.7	91.2	187.8
5	127.6126	0.3890	50.0	127.6	245.5	5	105.1525	0.4299	36.2	105.2	214.4
6	140.3556	0.3817	56.2	140.4	267.1	6	119.1360	0.4227	42.0	119.1	240.4
7	155.0333	0.3744	63.4	155.0	291.9	7	134.9677	0.4156	48.6	135.0	269.6
8	173.3780	0.3673	72.4	173.4	323.1	8	154.3755	0.4085	56.9	154.4	305.3
9	196.1267	0.3602	83.6	196.1	361.6	9	178.4732	0.4015	67.2	178.5	349.4
10	222.8130	0.3532	96.9	222.8	406.6	10	206.5732	0.3945	79.5	206.6	400.3
11	251.5695	0.3463	111.6	251.6	454.4	11	235.8698	0.3877	92.6	235.9	452.6
12	278.8827	0.3396	126.1	278.9	498.7	12	262.6764	0.3809	105.3	262.7	499.1
13	300.7811	0.3330	138.6	300.8	532.5	13	283.4969	0.3742	115.9	283.5	533.4
14	314.4302	0.3265	147.5	314.4	551.2	14	296.2184	0.3677	123.4	296.2	552.0

15	318.7692	0.3202	152.2	318.8	553.5	15	300.1497	0.3614	127.4	300.1	554.2
16	314.9775	0.3141	152.9	315.0	541.8	16	295.9452	0.3553	127.9	295.9	541.5
17	305.3809	0.3082	150.6	305.4	520.6	17	285.1932	0.3495	125.3	285.2	517.3
18	292.0343	0.3026	146.2	292.0	493.6	18	270.0078	0.3440	120.5	270.0	485.8
19	276.0573	0.2972	140.2	276.1	462.7	19	252.5661	0.3387	114.4	252.6	450.8
20	258.5452	0.2922	133.1	258.5	430.0	20	234.7593	0.3338	107.8	234.8	416.0
21-25	217.1600	0.2809	115.2	217.2	354.8	21-25	196.2189	0.3230	92.9	196.2	342.0
26-30	176.8353	0.2645	97.9	176.8	281.6	26-30	158.6612	0.3073	78.4	158.7	270.0
31-35	156.4434	0.2566	88.3	156.4	246.0	31-35	144.6606	0.2988	73.1	144.7	243.0
36-40	147.8521	0.2571	83.4	147.9	232.7	36-40	135.6504	0.2964	69.0	135.7	227.0
41-45	135.7125	0.2654	74.9	135.7	216.4	41-45	121.6626	0.2990	61.5	121.7	204.4
46-50	125.7431	0.2798	66.9	125.7	205.1	46-50	114.5344	0.3059	56.8	114.5	194.5
51-55	119.3813	0.2977	60.6	119.4	200.3	51-55	109.8845	0.3163	53.0	109.9	189.6
56-60	112.5445	0.3161	54.3	112.5	194.2	56-60	98.0102	0.3289	45.6	98.0	172.4
61-65	105.9676	0.3335	48.8	106.0	187.7	61-65	94.1497	0.3425	42.2	94.1	169.0
66-70	105.7704	0.3492	46.5	105.8	191.9	66-70	88.7193	0.3559	38.3	88.7	162.5
71-75	96.7081	0.3633	40.9	96.7	179.2	71-75	88.2148	0.3688	36.6	88.2	164.7
76-80	91.1025	0.3762	37.1	91.1	172.0	76-80	86.7241	0.3811	34.7	86.7	164.8
81-85	86.0783	0.3885	33.8	86.1	165.4	81-85	89.1368	0.3933	34.4	89.1	172.4
86-90	84.9649	0.4003	32.2	85.0	166.1	86-90	90.2972	0.4054	33.6	90.3	177.8

15. Performance Data

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1 Sensitivity

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 60 blanks and 100 low level samples.

Sensitivity	Concentration (ng/mL)
Limit of Blank (LoB)	1.9
Limit of Detection (LoD)	4.4
Limit of Quantitation (LoQ)	8.8

15.2 Precision

Precision was evaluated in accordance with a modified protocol based on CLSI EP-5A2, "Evaluation of Precision Performance of Quantitative Measurement Methods". A total of 5 serum samples were assayed using 3 lots of reagents in duplicate, twice a day for 20 days on one System.

Sample	N	Mean Conc. (ng/mL)	Within run		Total	
			SD	CV%	SD	CV%
1	80	21.9	0.6	2.9	1.2	5.4
2	80	81.2	1.9	2.3	2.6	3.5
3	80	163	3.1	1.9	6.4	3.9
4	80	663	15.5	2.3	35.0	5.3
5	80	304	4.2	1.4	21.8	7.2

15.3 Recovery

Recovery was assessed by adding known concentrations of IGF-I to low or medium level patient samples prior to assaying.

Sample Conc. (ng/mL)	Analyte added (ng/mL)	Measured (ng/mL)	Recovery (ng/mL)	Recovery (%)
92.0	120	221	129	108
92.0	240	356	264	110
88.0	120	208	120	100
88.0	240	328	240	100
178	120	297	119	99
178	240	419	241	100
Mean				103

15.4 Linearity

Linearity was evaluated based on CLSI EP-6A, "Evaluation of the Linearity of Quantitative Measurement Procedures". Samples were prepared by diluting a high patient sample with a low patient sample prior to assay. Each dilution is assayed in duplicate.

Predicted Concentration (ng/mL)	Measured Concentration (ng/mL)	Variation (ng/mL)	%
<10	11.9		
130	127	-3	-2%
256	266	10	4%
382	379	-3	-1%
509	475	-34	-7%
635	639	4	1%
761	779	18	2%
888	879	-9	-1%
1014	1022	8	1%

15.5 Method Comparison

The IDS-iSYS IGF-I assay was compared against a commercially available quantitative immunoassay, following CLSI EP-9A2, "Method Comparison and Bias Estimation Using Patient Samples". A total of 67 samples, selected to represent a wide range of IGF-I concentrations [15.5-613 ng/mL], was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

N	Slope	95% CI	Intercept (ng/mL)	95% CI	Correlation Coefficient (r)
67	0.62	0.57 to 0.67	-2.3	-4.5 to 9.1	0.99

15.6 Specificity

The specificity was assessed with the following cross-reactants.

Cross Reactant	% Cross-Reactivity
IGF-I	100 %
IGF- II	<0.01 %
Proinsulin	<0.01 %
Insulin	<0.01 %

16. Symbols used

Catalogue Number

In Vitro Diagnostic Device

Manufacturer

Applied in accordance with directive 98/79/EC

Caution: Federal law restricts this device to sale by or on the order of a physician (Only for US)

17. Bibliography

- Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. *Ann Clin Biochem* 1991; 28: 542-555.
- Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. *J Clin Endocrinol Metab* 1988; 67: 265-272.
- Barkan AL, Beitins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. *J Clin Endocrinol Metab* 1988; 67: 69-73.
- Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. *J Clin Endocrinol Metab* 1986; 67: 1348-1353.
- Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.
- Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. *J Clin Endocrinol Metab*. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idsplc.com • www.idsplc.com

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17. Bibliographie

1. Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. *Ann Clin Biochem* 1991; 28: 542-555.
2. Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. *J Clin Endocrinol Metab* 1988; 67: 265-272.
3. Barkan AL, Beittins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. *J Clin Endocrinol Metab* 1988; 67: 69-73.
4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. *J Clin Endocrinol Metab* 1986; 67: 1348-1353.
5. Boscart, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. *J Clin Endocrinol Metab*. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tél. : +44 191 519 0660 • Fax : +44 191 519 0760
courriel : info.uk@idsplc.com • www.idsplc.com

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17. Literatur

1. Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. *Ann Clin Biochem* 1991; 28: 542-555.
2. Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. *J Clin Endocrinol Metab* 1988; 67: 265-272.
3. Barkan AL, Beitins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. *J Clin Endocrinol Metab* 1988; 67: 69-73.
4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. *J Clin Endocrinol Metab* 1986; 67: 1348-1353.
5. Boscarto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. *J Clin Endocrinol Metab*. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
E-Mail: info.uk@idsplc.com • www.idsplc.com

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17. Bibliography

1. Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. Ann Clin Biochem 1991; 28: 542-555.
2. Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. J Clin Endocrinol Metab 1988; 67: 265-272.
3. Barkan AL, Beitzins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. J Clin Endocrinol Metab 1988; 67: 69-73.
4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. J Clin Endocrinol Metab 1986; 67: 1348-1353.
5. Boscaro LM and Stuart MC, 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. J Clin Endocrinol Metab. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idsplc.com • www.idsplc.com

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Toegepast conform richtlijn 98/79/EG

17. Bibliography

1. Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. Ann Clin Biochem 1991; 28: 542-555.
2. Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. J Clin Endocrinol Metab 1988; 67: 265-272.
3. Barkan AL, Beitins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. J Clin Endocrinol Metab 1988; 67: 69-73.
4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatotropin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. J Clin Endocrinol Metab 1986; 67: 1348-1353.
5. Boscaro LM and Stuart MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. J Clin Endocrinol Metab. 2014; 99(5): 1712-21.

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
E-mail: info.uk@idsplc.com • www.idsplc.com

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Výrobce

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17. Literatura

- Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. Ann Clin Biochem 1991; 28: 542-555.
- Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. J Clin Endocrinol Metab 1988; 67: 265-272.
- Barkan AL, Beitins IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. J Clin Endocrinol Metab 1988; 67: 69-73.

4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. *J Clin Endocrinol Metab* 1986; 67: 1348-1353.
5. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. *J Clin Endocrinol Metab*. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idspcl.com • www.idspcl.com

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17. Bibliography

1. Rudd BT. Growth, growth hormone, and the somatomedins: A historical perspective and current concepts. *Ann Clin Biochem* 1991; 28: 542-555.
2. Baxter RC. Characterisation of the acid-labile subunit of the growth hormone dependant insulin-like growth factor binding protein complex. *J Clin Endocrinol Metab* 1988; 67: 265-272.
3. Barkan AL, Beutels IZ, Kelch RP. Plasma immunoreactive insulin-like growth factor 1 in acromegaly: Correlation with the degree of growth hormone hypersecretion. *J Clin Endocrinol Metab* 1988; 67: 69-73.
4. Oppizi G, Petroncini M, Dallabonzana D, et al. Relationship between somatomedin C and growth hormone levels in acromegaly : Basal and dynamic evaluation. *J Clin Endocrinol Metab* 1986; 67: 1348-1353.
5. Boscarto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.
6. Bidlingmaier M, Friedrich N, Emeny RT, et al. Reference Intervals for Insulin-like Growth Factor-1 (IGF-1) From Birth to Senescence: Results From a Multicenter Study Using a New Automated Chemiluminescence IGF-1 Immunoassay Conforming to Recent International Recommendations. *J Clin Endocrinol Metab*. 2014; 99(5): 1712-21.

 Immunodiagnostic Systems Limited,

10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idsplc.com • www.idsplc.com

Immunodiagnostic Systems

EUA Immunodiagnostic Systems (IDS) Inc,
948 Clopper Road, Gaithersburg, MD 20878, USA.
Tel.: +1 877 852 6210 • Fax: +1 301 990 4236
e-mail: info.us@idsplc.com

França Immunodiagnostic Systems SAS,
153 Avenue D'Italie, 75013 Paris, France
Tel.: +33 140 770 450 • Fax: +33 140 770 455
e-mail: info.fr@idsplc.com

Bélgica Immunodiagnostic Systems S.A.,
Rue Ernest Solvay 101, 4000 Liège, Belgium
Tel.: +32 425 226 36 • Fax: +32 425 251 96
e-mail: info.be@idsplc.com

Alemanha Immunodiagnostic Systems GmbH (IDS GmbH),
Rahmhofstraße 2-4, 60313 Frankfurt am Main
Tel.: +49 69 26019 09650 • Fax: +49 69 26019 0949
e-mail: info.de@idsplc.com

Escandinávia Immunodiagnostic Systems Nordic a/s (IDS Nordic a/s), International House, Center Boulevard 5,
2300 København S, Denmark
Tel.: +45 448 400 91
e-mail: info.nordic@idsplc.com

Brasil IDS Brasil Diagnósticos Ltda.
Rua dos Pinheiros, 610 - conj 41 Edifício Win Work Pinheiros
São Paulo, SP Brasil. CEP: 05422 001
Tel.: +55 113 740 6100 • Fax: + 55 113 740 6105
e-mail: info.br@idsplc.com