SUHRN A VYSVETLENIE

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SUMMARY AND EXPLANATION

VP reagent is auxiliary reagent for test VP (Voges – Proskauer test). Reagent shows colour change by creating of acetoín in tests of diagnostic kits and strips produced by DIAGNOSTICS s.r.o.

PRINCIPLE

Glucose is fermented by microorganisms to acetoín. Reagent VP, 1 – ethanol solution of α -naftol, deepens final colouring of solution and simultaneously elevates accuracy of reaction. α - naftol catalyses change of acetoín in alkaline environment (VP2) to diacetyl. KOH reacts with peptone. Condensation of diacetyl with guanid groups of peptone in alkaline environment leads to formation of pink colour, what represents positive result of test.

DIAGNOSTIC KIT CONTAINS

- reagent VP 1a – substrate (Ref. 3004)
- reagent VP 1b – solvent (Ref. 3004)
- reagent VP 2 (Ref. 3004)
- dropper
- information leaflet

PREPARING OF REAGENT TO USE

- pour VT 1a to VP1b. After dissolving of substrate, reagent is ready to use. Mark the label of reagent with the preparation and expiration date. Reagent VP 2 is ready for use.

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety regulations of your country must be respected while handling
- Do not use overexpired product.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

Kit contains ethanol (50 % of volume), methanol (3% of volume), hydroxid draselný (20 % of volume) , α – naftol (< 1 % of volume)

R11: Highly flammable.
R20: Harmful by inhalation.
R35: Causes burns
R 36/R37/R38: Irritating eyes, respiratory system and skin.
S 26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice.
S 36/37/39 Wear suitable protective equipment.
S 45 In case of accident or if you feel unwell seek medical advice immediately (show the label when possible).

STORAGE

Store VP reagent in dark, in temperatures between +2 and +8°C. Expiration date is showed on label of each package. Prepared reagent VP1 (VP 1a to VP 1b) store in dark at temperatures between +2 and +8°C. Use prepared reagent till 5 months from date of preparation, or till date of expiration, - in case, that expiration date pass earlier.

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