diagnostics HIP reagent

Pomocné činidlo / Auxiliary reagent / Pomocné činidlo / Segédreagens

IVD

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Dátum poslednej revízie: 8.1.2014

SUMMARY AND EXPLANATION

HIP reagent is reagent used for generation of colour in test HIP in diagnostic kit. Detailed description of work with reagent HIP and interpretation of result of test is part of particular kit manual.

PRINCIPLE

Substrate sodium hippurate is hydrolysed to end product, glycine and benzoic acid, which is detected by HIP reagent.

DIAGNOSTIC KIT CONTAINS

- HIP solution 20 ml (a) (Ref. 3006)
- HIP reagent - solid for reagent (b) (Ref. 3006)
- dropper
- information leaflet

PREPARING OF REAGENT TO USE

Pour solution "a" to solid "b", then put on dropper. After dissolving of solid, reagent is prepared to be used.

WARNINGS AND SECURITY PRECAUTIONS

- For *in vitro* diagnostics use.
- For professional use only
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of overexpired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

Kit contains mixture of acetone/butanol (50 % / 50 %) and ninhydrin (100 % of solid).

R11: Highly flammable.

R20: Harmful by inhalation.

R22: Harmful if swallowed.

R 36/R37/R38: Irritating eyes, respiratory system and skin.

S 26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice.

S 36 Wear suitable protective clothing.

S 45 In case of accident or if you feel unwell seek medical advice immediately (show the label when possible).

STORAGE

HIP reagent can be stored in dark at temperatures between +2 and +8°C. Expiration date is showed on label of each package.

Reagent prepared by mixing of components can be stored in dark in temperatures between + 2°C and + 8°C for 5 months or till date of expiration – if expiration date finishes earlier.

Warning: When preparing solution of ninhydrin avoid contact with skin. Contact leads to formation of spots that are hard to remove.

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