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Dátum poslednej revízie: 7.5.2014

SUMMARY AND EXPLANATION

Diagnostic kit CARBA is intended for fast identification of bacteria producing carbapenemases.

PRINCIPLE

By enzymatic hydrolysis of β -lactam ring, colouring of inoculated solution changes from negative red-pink to yellow, indicating positive test result of examined microorganism.

DIAGNOSTIC KIT CONTAINS

- CARBA „a“ - substrate in 5 pc of microtubes eppendorf
- CARBA „b“ - solution 5 ml
- information leaflet

RECOMMENDED PROCEDURE

- Add 12 - 14 drops of solution CARBA „b“ to microtube CARBA „a“. Put one full 10 μ l inoculation loop of bacteria mass to prepared microtube and thoroughly homogenise. Use fresh 24 h bacterial culture from suitable cultivation medium (for example blood agar).
- Evaluate result of test after 0,5-2 hours.

For the need of own proof function, use following recommended bacterial strains:

Bacterial strain	CCM / ATCC	Colour	Interpretation
<i>Salmonella enterica subsp. enterica</i>	4419 / -	yellow	positive
without bacteria	-	red-pink	negative

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

WARNINGS AND SECURITY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Store CARBA at dark place at +2 to +8 °C. Expiration date is placed on the label of each package.

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