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Dátum poslednej revízie: 10. 3. 2014

SUMMARY AND EXPLANATION

PYR is intended for the detection of pyrrolidonylarylamidase enzymatic activity. Detection strip performs the colour change by creating of orange or red colour based on the metabolism of investigated samples. The strip is used for confirmation of *Streptococcus pyogenes* and species belonging to the genus *Enterococcus*.

PRINCIPLE

The substrate is broken down by enzymatic activity of investigated microorganism, what result in formation of red colour on the surface of diagnostic zone after adding PYR reagent (Ref. 3003) (containing p-dimethylaminocinnamaldehyde).

DIAGNOSTIC KIT CONTAINS

- 50 test strips for 50 determinations
- information leaflet

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious valid regulations must be respected when handling.
- Observe common safety measures according to the regulations applicable in each country.
- Before use, check the packaging is intact. Do not use damaged kit.

STORAGE

Store in the dark at +2 to +8°C. Expiration date is placed on the label of each package.

REQUIRED MATERIALS

PYR reagent (Ref. 3003)

RECOMMENDED PROCEDURE

A) spot test option

1. Use pure 24 hour bacterial culture. Smear several well separated colonies on diagnostic zone of the strip. If necessary, moisten the zone with approximately 10 µl of water.
2. Incubate at the appropriate plate for 5-10 minutes.
3. Add drop of PYR reagent after incubation and read the reaction visually after 1-2 minutes.

B) tube test option

1. Use pure 24 hour bacterial culture. Take a tube with 0,5 ml saline. Prepare a homogenous suspension with density of turbidity for at least 3 McF.
2. Insert the detection strip and dip detection zone into bacterial suspension and incubate at 35±2°C for 2-4 hours.
3. Add 3 drops of PYR reagent into the tube after incubation and evaluate result of test.

INTERPRETATION OF REACTION

Red or orange colour indicates positive result of the PYR test.

For the need of own proof function, use following recommend bacterial strains:

Bacterial strain	CCM / ATCC	PYR	Colour
<i>Enterococcus faecalis</i>	4224 / 29212	+	red, orange

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

TIP: use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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