VP

Detekčný prúžok pre mikrobiológiu / Diagnostic strip for microbiology / Detekční proužek pro mikrobiologii / In vitro diagnostikum

IVD

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Dátum poslednej revízie: 10. 3. 2014

SUMMARY AND EXPLANATION

VP is intended for the detection of acetoine production in Voges-Proskauer reaction. Detection strip performs the colour change creating pink or red colour based on the metabolism of investigated samples.

PRINCIPLE:

Glucose is fermented by some microorganisms to acetoine. VP 1 reagent, ethanol solution of α -naphthol, deepens the final colouring and increases the sensitivity of the reaction. α -naphthol catalyses the conversion of acetoine in an alkaline medium (VP 2) to form diacetyl. KOH reacts with peptone. Diacetyl condenses with guanidin groups of peptone leading to formation of pink or red colour, what indicates positive result of the test.

DIAGNOSTIC KIT CONTAINS:

- 50 test strips for 50 determinations
- information leaflet

WARNINGS AND SAFETY PRECAUTIONS

- For *in vitro* diagnostics use.
- For professional use only.
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious valid regulation must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

STORAGE

In refrigerator at +2 to +8°C. Store in the dark. Expiration date is placed on the label of each package.

REQUIRED MATERIALS

VP reagent (Ref. 3004)

Tube with 0,5-1 ml saline

RECOMMENDED PROCEDURE:

1. Use pure 24 hour bacterial culture. Take a tube with 0,5 ml saline. Prepare a homogenous suspension with density turbidity from 2 to 3 McF.
2. Insert the detection strip and dip detection zone into bacterial suspension.
3. Incubate at 35±2°C for 2 hours.
4. Add 3 drops of VP1 reagent and then 3 drops of VP2 reagent into the tube after incubation.
5. Incubate at 35±2°C for 30 minutes and evaluate result of a test.

INTERPRETATION OF REACTION

Red or pink colour indicates positive result of the VP test.

For the need of own proof function, use following recommend bacterial strains:

Bacterial strain	CCM / ATCC	Test	Colour
<i>Staphylococcus aureus</i>	4223 / -	positive	red, pink

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420 549491430, e-mail: ccm@sci.muni.cz

TIP: use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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