

FUCAST-and CLSI potency NEO-SENSITABS™

Page 1 / 4 - Document: 3.15.0

Interpretation Zones and MIC Breakpoints according to CLSI

Yeasts (CLSI's M44-A2 Method of Diffusion on Agar)

Only Fluconazole and Voriconazole are included in the M44-A guideline from CLSI (formerly NCCLS)(1) and interpretation zones are only available for Fluconazole so far. When MIC breakpoints are available from CLSI (2) they are used in the table below.

Table 13.5-1 Interpretation for Yeasts

Mueller-Hinton Glucose Methylene Blue agar. Inoculum McFarland 0.5 undiluted. Incubation at 35 °C ± 2 °C ambient air for 20-24 hours. MICs according to M27-S3 (2007).

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Amphotericin B	10 µg	AMPHO	≥ 15	14-10	< 10	≤ 1	≥ 2
Caspofungin *) (24)	5 µg	CASP5					
<i>C. albicans</i> , <i>C. tropicalis</i>			≥ 17	16-15	≤ 14	≤ 0.25	≥ 1
<i>C. parapsilosis</i> ,* <i>C. guilliermondii</i>			≥ 13	12-11	≤ 10	≤ 2	> 2
<i>C. glabrata</i> , <i>C. krusei</i>			≥ 16	-	-	≤ 0.5	-
Fluconazole **) (27)	25 µg	FLUCZ					
<i>C. albicans</i> , <i>C. tropicalis</i> , <i>C. parapsilosis</i>			≥ 17	16-14	≤ 13	≤ 2	≥ 8
<i>C. glabrata</i>			-	≥ 15	≤ 14	≤ 32	≥ 64
(SDD)						(SDD)	
Itraconazole (11)	10 µg	ITRAC	≥ 23	22-14	≤ 13	≤ 0.12	> 0.5
(SDD)							
Ketoconazole	15 µg	KETOC	≥ 28	27-21	≤ 20	≤ 0.12	≥ 0.5
Posaconazole (9,12)	5 µg	POSAC	≥ 17	16-14	≤ 13	≤ 1	> 2
(SDD)							
Voriconazole (23)	1 µg	VOR.1					
<i>C. albicans</i> , <i>C. tropicalis</i> , <i>C. parapsilosis</i>			≥ 17	16-15	≤ 14	≤ 0.12	≥ 1
<i>C. krusei</i>			≥ 15	14-13	≤ 12	≤ 0.5	≥ 2
<i>C. glabrata</i>			≥ 16	-	≤ 15	≤ 0.5	> 0.5

*) Tentative. There is cross-resistance between Caspofungin and the other echinocandins: Anidulafungin and Micafungin (19).

**) *C. krusei* should be reported as resistant to Fluconazole (no matter the zone).

For further information on Susceptibility Testing of Yeasts, see Neo-Sensitabs User's Guide.

Table 13.5-2 Mold Disk Diffusion Testing (14, 17,21)

Mueller-Hinton Plain. Incubation for 16-24 hours (zygomycetes), 24 hours (Aspergillus), 48 hours (other spp.). Temperature: 35-37°C

Inoculum:

Prepare a suspension of sporulating colonies in 0.85% saline, add 1 drop of Tween 20. Allow heavy particles to settle for 3-5 minutes and the upper suspension is treated for 15 seconds in a vortex mixer. The density of the suspension is read on a spectrophotometer at 530 mm wave length and the optical density adjusted at 0.09 to 0.13 for Aspergillus.

**Interpretation Zones and MIC
Breakpoints according to CLSI**

**Yeasts (CLSI's M44-A2 Method of Diffusion
on Agar)**

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Amphotericin B (zygomycetes only)	10 µg	AMPHO	≥ 15	-	-	≤ 1	-
Caspofungin	5 µg	CASP5	≥ 17	-	-	≤ 1	-
Itraconazole	10 µg	ITRAC	≥ 17	-	-	≤ 1	≥ 2
Posaconazole	5 µg	POSAC	≥ 23	-	-	≤ 0.25	≥ 0.5
Voriconazole	1 µg	VOR.1	≥ 17	-	-	≤ 1	≥ 2

Note: The base medium should not be supplemented with neither 2% glucose nor 0.5% methylene blue dye.

Itraconazole susceptible strains should be reported as susceptible to both, Posaconazole and Voriconazole.

Interpretation table for Local treatment

In local treatment of fungal infections, a high concentration of antifungal is placed at site of the infection. Consequently other MIC breakpoints and zone interpretations should be used in those cases.

Local Treatment MH Glucose Methylen Blue Agar or Shadowy McFarland 0.5 inoculum			
Susceptible	≥ 20 mm	≥ 15 mm	≥ 10 mm
Intermediate	12-19 mm	10-14 mm	-
Resistant	≤ 11 mm	no zone	no zone
	Ciclopirox Clotrimazole Econazole Fluconazole Isoconazole Ketoconazole Miconazole Tioconazole Terbinafine	Natamycin Nystatin Itraconazole	Griseofulvin

Fluorocytosine cannot be tested on MH-agar (antagonists), but has to be tested on Shadowy agar or similar.

Table 13.5-3 Quality Control Zone Diameters (mm) Ranges

Mueller-Hinton Glucose Methylene Blue agar. Inoculum McFarland 0.5 undiluted. Incubation at 35 °C ± 2 °C for 20-24 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm		
			<i>C. albicans</i> ATCC 90028	<i>C. parapsilosis</i> ATCC 22019	<i>C. krusei</i> ATCC 6258
Amphotericin B	10 µg	AMPHO	20-27	22-29	18-25
Fluconazole	25 µg	FLUCZ	28-39	22-33	-
Itraconazole	10 µg	ITRAC	21-30	19-26	16-22
Ketoconazole	15 µg	KETOC	31-42	35-45	22-29
Voriconazole	1 µg	VOR.1	31-42	28-37	23-31
Caspofungin	5 µg	CASP5	15-22	13-23	16-22
Posaconazole	5 µg	POSAC	24-34	25-36	23-31

Carrillo-Munoz et al (30) determined Posaconazole susceptibility of clinical yeast isolates with Neo-Sensitabs and a microdilution method. Complete agreement between Posaconazole Neo-Sensitabs and the microdilution was 92.3 % after 24 hours incubation. The authors conclude that the Agar Diffusion with Posaconazole Neo-sensitabs can improve Posaconazole susceptibility testing due to its excellent correlation and reduced percentage of disagreements in comparison with microdilution testing.

References:

- 1) NCCLS. Methods for Antifungal Disk Diffusion Susceptibility Testing on Yeasts. Approved Guideline M44-A2, 2008.
- 2) NCCLS. Reference Methods for Broth Dilution Antifungal Susceptibility Testing of yeasts. Approved Standard, 2nd Ed. M27-A2, 2002.
- 3) Vandebossche I. et al: Susceptibility testing of Fluconazole by the NCCLS broth macrodilution method. E-Test and Disk diffusion for application in the routine laboratory. J. Clin. Microbiol., **40**, 918-921, 2002.
- 4) Pfaller M.A. et al: Evaluation of the NCCLS M44-P Disk diffusion method for determining Fluconazole susceptibility of 276 clinical isolates of *Cryptococcus neoformans*. ICAAC 2003, presentation M-1204.
- 5) Barry A. et al: Quality Control for Fluconazole disk susceptibility tests on Mueller-Hinton Agar with glucose and methylene blue. J.Clin.Microbiol., **41**, 3410-12, 2003.
- 6) Rementeria A. et al: Utility of Neo-S. tablets of Fluconazole and Voriconazole for in vitro susceptibility testing of *Candida* spp. with the NCCLS M44-P method of diffusion on agar. VII Congress Mycology. Salamanca, July 2004 (Spanish).
- 7) Espinel-Ingroff A. et al: Correlation between Neo-S. tablets on 3 media, NCCLS reference disk diffusion and broth microdilution methods for testing *Candida* spp. and *Cryptococcus neoformans* with Fluconazole and Voriconazole. Abstract, ECCMID 2005.
- 8) Espinel-Ingroff A, Canton E.: Comparison of 3 antifungal susceptibility methods for testing *Candida* spp. & *Cryptococcus neoformans* with Caspofungin (CAS) and Amphotericin B (AMB). Abstract, presentation M-1604, 45th ICAAC 2005.
- 9) Sims C.R. et al: Correlation between microdilution, E-test and disk diffusion methods for antifungal susceptibility testing of Posaconazole against Candida spp. J. Clin. Microbiol., **44**, 2105-08, 2006.
- 10) Brown S. et al: Q.C. limits for Posaconazole disk susceptibility tests on MH agar with glucose and methylene blue. J. Clin. Microbiol., **45**, 222-23, 2007.
- 11) Espinel-Ingroff A. et al: Correlation of Neo-Sensitabs tablet diffusion assay results on 3 different agar media with CLSI broth microdilution M27-A2 and disk diffusion M44-A results for testing

- susceptibilities of *Candida* spp. and *Cryptococcus neoformans* to Amphotericin B, Caspofungin, Fluconazole, Itraconazole and Voriconazole. *J. Clin. Microbiol.*, **45**, 858-64, 2007.
- 12) Diekema D.J. et al: Evaluation of E-test and disk diffusion methods with broth microdilution antifungal susceptibility testing of Posaconazole against clinical isolates of *Candida*. *J. Clin. Microbiol.*, 2007 (in press).
- 13) Carrillo-Muñoz A.J. et al: Acitivity of Caspofungin and Voriconazole against clinical isolates of *Candida* and other medically important yeasts by the CLSI M44A disk diffusion method with Neo-Sensitabs tablets. *Chemotherapy*, **54** (1) 38-42, 2008.
- 14) Espinel-Ingroff A. et al: Multicenter evaluation of a new disk agar diffusion method for susceptibility testing of filamentous fungi with Voriconazole, Posaconazole, Itraconazole, Amphotericin B and Caspofungin. *J. Clin. Microbiol.*, **45**, 1811-20, 2007.
- 15) Espinel-Ingroff A.: Standardized disk diffusion method for yeasts. *Clin. Microbiol. Newsl.*, **29**, 97-100, 2007.
- 16) Carrillo-Muñoz A. J. et al: Antifungal activity in vitro of Posaconazole. Method of diffusion on agar (Neo-Sensitabs) based on document M44-A, against yeast. III. Forum Micologico, Dec 15th 2007, Alicante, Spain (spanish).
- 17) Espinel Ingroff A. et al: Comparison of Neo-S tablet diffusion assay with CLSI broth microdilution M38-A and disk diffusion methods for testing susceptibility filamentous fungi with Amphotericin B, Caspofungin, Itraconazole, Posaconazole and Voriconazole. *J. Clin. Microbial.* **46**, 1793-1803, 2008.
- 18) Duran-Valle M T et al: Resistance to Voriconazole and/or Fluconazole in yeast isolated in a Madrid hospital (2005-2007) etc. XII Congress SEIMC, presentation 33, Madrid, May 2008.
- 19) Pfaller M. A. et al: Correlation of MIC with outcome for *Candida* spp. tested against Caspofungin, Amidulafungin and Micafungin: analysis and proposal for interpretative MIC breakpoints. *JCM* **46**, 2620-2629, 2008.
- 20) Pfaller M. A. et al: Wild-type MIC distribution and epidemiological cutoff values for *Aspergillus fumigatus* and 3 triazoles as determined by the CLSI both microdilution methods. *J. Clin. Microbiol.* **47**, 3142-46, 2009.
- 21) M51-A Method for antifungal disk diffusion susceptibility testing of filamentous fungi approved guideline. CLSI 2009.
- 22) Pfaller M. et al: Clinical breakpoints for the Echinocandins and *Candida*, revisited. Presentation M-369, 50th ICAAC, Sept. 2010.
- 23) Pfaller M. et al: Clinical breakpoints for Voriconazole and *Candida* spp. revisited: Review of microbiologic, molecular, pharmacodynamic and clinical data as they pertain to the development of species-specific interpretative criteria. *Diagn. Microbiol. Infect. Dis.* **70**, 330-343, 2011.
- 24) Pfaller M. et al: Clinical breakpoints for Echinocandins and *Candida* revisited: Integration of molecular, clinical and microbiological data to arrive at species-specific interpretative criteria. *Drug Resistance Updates*. **14**, 164-176, 2011.
- 25) Ochiuzzi M. E et al: Correlation of E-test and Neo-Sensitabs diffusion assays on Mueller-Hinton-methylene blue agar, with broth microdilution reference method (CLSI-M27-A2) for testing susceptibilities of *Cryptococcus neoformans* to Amphotericin B and Fluconazole. *Med. Mycology*. **48**, 893-896, 2010.
- 26) Carrillo-Muñoz A. J. et al: Antifungal activity of Posaconazole against *Candida* spp. and non-*Candida* clinical yeast isolates. *Rev. Esp. Quimioter.* **23**, 122-125, 2010.
- 27) Nweze E. I et al: Agar-based disk diffusion assay for susceptibility testing of Dermatophytes. *J. Clin. Microbiol.* **48**, 3750-52, 2010.
- 28) Pfaller M. A et al: Wild-type MIC distributions, epidemiological cutoff values and species-specific clinical breakpoints for Fluconazole and *Candida*: Time for harmonization of CLSI and EUCAST broth microdilution methods. *Drug Resistance Updates*. **13**, 180-195, 2010.
- 29) Pfaller M. A et al: Echinocandin and triazole antifungal susceptibility profiles for *Candida* spp., *Cryptococcus neoformans* and *Asp. fumigatus*; application of new CLSI clinical breakpoints and epidemiologic cutoff values to characterize resistance in SENTRY Antimicrobial Surveillance Program (2009). *Diag. Microbiol. Infect. Dis.* **69**, 45-50, 2011.
- 30) Carrillo-Munoz AJ et al: Posaconazole susceptibility of clinical yeast isolates determined by an agar diffusion and microdilution method. *Int. J. Antimicr. Agents*. **37**, 271-273, 2011.