

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test Serum Anti-Coli O 26:K 60**

Catalogue No. (Package size) TS 2202 (1 ml)

Lot No. **104 07 13**

Expiry date **2014-03**

Intended use Test Serum Anti-Coli O26:K60 (undiluted) records all E. coli strains that contain the antigen O26:K60 in slide agglutination. In the Widal reaction the serum reaches an O titre of 1:800 (tube) or 1:1600 (MTP) with the homologous E. coli strains; deviations of ± 1 titre step are acceptable.

Specificity Since unspecific co-agglutinations with other E. coli strains often occur in the slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled E. coli isolate in the Widal test.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin, 23.07.13


Head of Quality Control