

np nerbe plus

GERMANY

SEROLOGISCHE PIPETTEN / SEROLOGICAL PIPETTES



ICONS

ICONS



Artikelnummer
Order no.



Information
Information



Material
Material



Neues Produkt
New product



Innenkarton
Box



Abmessungen
Measurements



Umkarton
Master shipper



Steril R (strahlensteril)
Sterile R (sterile by radiation)



Polybeutel
Bag




mp nerbe plus

10 ml in 1/10
STERILE IF UNOPENED
AND UNDEFORMED

10ml in 1/10
nerbe plus
10ml in 1/10
nerbe plus
10ml in 1/10
nerbe plus

SEROLOGISCHE PIPETTEN SEROLOGICAL PIPETTES

- **Fussel freier Filter aus PE**
Lint-free PE filter
- **Nicht zytotoxisch**
Non-cytotoxic
- **Farbcodiert**
Color coded



12-521-9101

12-321-6103

12-331-6103

12-401-9103

12-401-6103

12-421-9101

12-421-6101

REF	MAT						
12-521-9101	PS	275 mm	2 ml	50		500	Aspirationspipette, Einzel-Peelpack, ohne Farbcode / aspirating pipette, single peelpack, no color code STERILE R
12-321-6103	PS	270 mm	1,1 ml	25		1.000	Farbcode: orange / color code: orange STERILE R
12-331-6103	PS	275 mm	2,2 ml	25		1.000	Farbcode: grün / color code: green STERILE R
12-401-6103	PS	270 mm	1 ml	25		1.000	Farbcode: gelb / color code: yellow STERILE R
12-401-9103	PS	270 mm	1 ml			500	Einzel-Peelpack, Farbcode: gelb / single peelpack, color code: yellow, Dispenser-Box STERILE R
12-421-6101	PS	275 mm	2 ml	25		1.000	Farbcode: grün / color code: green STERILE R
12-421-9101	PS	275 mm	2 ml			500	Einzel-Peelpack, Farbcode: grün / single peelpack, color code: green, Dispenser-Box STERILE R



np sterile plus
Serologische Pipetten
Serological pipettes

Produktionsnr.: 81100-0001-0001
Nicht sterilisierbar / Not sterilizable
Nicht reusen / Do not reuse

STERILE IF UNOPENED OR
UNDAMAGED PACKAGE

Modell	Volumen	Skala	Skalierung	Skalierung	Skalierung
np 1000-1000	1000 µl	0,1	0,2	0,5	1
np 2000-1000	2000 µl	0,1	0,2	0,5	1
np 5000-1000	5000 µl	0,1	0,2	0,5	1
np 10000-1000	10000 µl	0,1	0,2	0,5	1

2 ml

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus

np sterile plus


np sterile plus

np sterile plus

SEROLOGISCHE PIPETTEN

SEROLOGICAL PIPETTES

- **Fusselfreier Filter aus PE**
Lint-free PE filter
- **Nicht zytotoxisch**
Non-cytotoxic
- **Farbcodiert**
Color coded



12-441-6105
12-441-9105
12-441-9115

12-461-6108
12-461-9108
12-461-9118


12-481-9102
12-481-9112

12-491-9102

REF	MAT						
12-441-6105	PS	335 mm	5 ml	25		500	Farbcode: blau / color code: blue STERILE R
12-441-9105	PS	335 mm	5 ml			200	Einzel-Peelpack, Farbcode: blau / single peelpack, color code: blue, Dispenser-Box STERILE R
+ 12-441-9115	PS	335 mm	5 ml	25		300	Einzel-Peelpack, Farbcode: blau / single peelpack, color code: blue STERILE R
12-461-6108	PS	340 mm	10 ml	25		500	Farbcode: orange / color code: orange STERILE R
12-461-9108	PS	340 mm	10 ml			200	Einzel-Peelpack, Farbcode: orange / single peelpack, color code: orange STERILE R
+ 12-461-9118	PS	340 mm	10 ml	25		300	Einzel-Peelpack, Farbcode: orange / single peelpack, color code: orange, Dispenser-Box STERILE R
12-481-9102	PS	345 mm	25 ml			150	Einzel-Peelpack, Farbcode: rot / single peelpack, color code: red, Dispenser-Box STERILE R
+ 12-481-9112	PS	345 mm	25 ml	25		200	Einzel-Peelpack, Farbcode: rot / single peelpack, color code: red STERILE R
12-491-9102	PS	345 mm	50 ml			100	Einzel-Peelpack, Farbcode: lila / single peelpack, color code: violet, Dispenser-Box STERILE R

SEROLOGISCHE PIPETTEN SPEZIALMODELLE

SEROLOGICAL PIPETTES SPECIAL



REF	MAT						
12-441-7105	PS	230 mm	5 ml	50		200	Einzel-Peelpack, kurze Ausführung , Farbcode: blau / single peelpack, short version , color code: blue STERILE R
12-461-7108	PS	235 mm	10 ml	25		200	Einzel-Peelpack, kurze Ausführung , Farbcode: orange / single peelpack, short version , color code: orange STERILE R
12-471-8108	PS	340 mm	10 ml	25		500	große Öffnung , Farbcode: orange / big orifice , color code: orange STERILE R
12-471-9108	PS	340 mm	10 ml			200	Einzel-Peelpack, große Öffnung , Farbcode: orange / single peelpack, big orifice , color code: orange STERILE R
12-481-7105	PS	230 mm	25 ml	25		100	Einzel-Peelpack, kurze Ausführung , Farbcode: rot / single peelpack, short version , color code: red STERILE R



nerbe plus GmbH & Co. KG
Porschestraße 25
21423 Winsen/Luhe · Germany
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 645 45
E-Mail: info@nerbe-plus.de
Web: www.nerbe-plus.de

