

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

|                    |   |
|--------------------|---|
| <b>Product</b>     | <b>Anti-Coli O 103:K-</b>   |
| Catalogue no.      | TS2216  |
| Package size       | 1 ml  |
| <b>Lot no.</b>     | <b>8961123</b>  |
| <b>Expiry date</b> | <b>2028-10</b>  |
| Intended use       | Anti-Coli O 103:K- records all <i>E. coli</i> strains that contain the antigen O 103:K- in slide agglutination. In the Widal reaction the serum reaches a titre of 1:3200 (tube) and 1:6400 ( MTP) with the homologous <i>E. coli</i> strains; deviations of $\pm 1$ titre step are acceptable. |
| Specificity        | Since unspecific co-agglutinations with other <i>E. coli</i> strains often occur in the slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled <i>E. coli</i> isolate in the Widal test.   |
| Preservatives      | Sodium azide 0.9 mg/ ml   |
| Storage            | 2 °C...8 °C   |

**This document has been established electronically and is valid without signature.  
The batch release was done in accordance with SOP QMV 7.**