

# XSX 200

Working Volume	With Screw Caps	220 µl
	With Push Caps	220 µl
	With Heat Seal	250 µl
Type of Rack	Polycarbonat	
Type of Rack Lid	Slide Lock	
Rack Identifier	Narrow Side	1D and Human Readable
	Bottom	2D
Storage Temperature	With Screw Caps	-196°C
	With Push Caps	-80°C
	With Heat Seal	-80°C
Standard Side Code Tube	N.A.	
Screw Caps	Pre-Capped	
	Stacked on SBS 96 Cap Rack Bulk Ware	
Purity	Production Sterile (ISO 8) Pre Sterilized (Beta Radiation)	
Outer Diameter	8.65 mm	
Tube Height	Without Cap	14.6 mm
	With Screw Cap	16.0 mm
	With Push Cap	16.0 mm
Rack Height	Without Cap	15.1 mm
	Screw Capped Tubes	16.5 mm
	With Lid	22.3 mm



## Technical Drawing



Standard Item Number 2DNC-X02-BL-NS-SLC-L

Order Information

