


SX 300

Working Volume	With Screw Caps	300 μ l
	With Push Caps	300 μ l
	With Heat Seal	330 μ l
Type of Rack	Polycarbonate Polypropylene	
Type of Rack Lid	Slide Lock	
Rack Identifier	Narrow Side	1D and Human Readable
	Bottom	2D
Storage Temperature	With Screw Caps	-196°C
	With Push Caps	-80°C
	With Heat Seal	-80°C
Standard Side Code Tube	(optional)	Human Readable
Screw Caps	Pre-Capped Stacked on SBS 96 Cap Rack Bulk Ware	
Purity	Production Sterile (ISO 8) Pre-Sterilized (Beta Radiation)	
Outer Diameter	8.65 mm	
Tube Height	Without Cap	18.2 mm
	With Screw Cap	19.5 mm
	With Push Cap	19.5 mm
Rack Height	Without Cap	18.7 mm
	Screw Capped Tubes	20.2 mm
	With Lid	22.3 mm



Technical Drawing



Standard Item Number 2DNC-X03-BL-NS-SLC-L

Order Information

