



MX 500

Working Volume	With Screw Caps	500 µl
	With Push Caps	500 µl
	With Heat Seal	530 µl
Type of Rack	Polycarbonate Polypropylene	
Type of Rack Lid	Slide Lock Push Lock	
Rack Identifier	Narrow Side	1D and Human Readable
	Bottom	2D
Storage Temperature	With Screw Caps	-196°C
	With Push Caps	-80°C
	With Heat Seal	-80°C
Standard Side Code Tube	(optional)	1D and Human Readable
Screw Caps	Pre-Capped Stacked on SBS 96 Cap Rack Bulk Ware	
Purity	Production Sterile (ISO 8) Pre-Sterilized (Beta Radiation)	
Outer Diameter	8.65 mm	
Tube Height	Without Cap	26.2 mm
	With Screw Cap	27.5 mm
	With Push Cap	27.5 mm
Rack Height	Without Cap	26.7 mm
	Screw Capped Tubes	28.0 mm
	With Lid	30.3 mm



Technical Drawing



Standard Item Number 2DNC-X05-BL-NS-SLC-L

Order Information

