

## Certificate of Quality

**Product Name:** VERSATREK MYCO GS  
**Lot Number:** 976700

**Product Number:** 7112-42  
**Expiration Date:** 2022-05-04  
(YYYY-MM-DD)

This product has been manufactured, processed and packaged in accordance with Quality Systems Regulation, 21 CFR Part 820. Representative samples were tested per Remel Inc., a part of Thermo Fisher Scientific Quality Control specifications and were found to meet performance criteria for this product.

<b>Organism Challenge</b>	<b>ATCC</b>	<b>Performance Criteria</b>
Mycobacterium fortuitum	6841	Positive in $\leq$ 14 Days
Mycobacterium intracellulare	13950	Positive in $\leq$ 14 Days
Mycobacterium kansasii	12478	Positive in $\leq$ 14 Days
Mycobacterium scrofulaceum	19981	Positive in $\leq$ 14 Days
Mycobacterium tuberculosis	25177	Positive in $\leq$ 14 Days

Product was tested in the Versatrek system as documented in the instructions for use.

### **Physical Characteristics:**

Appearance: Light Amber, Clear Broth w/ No Precipitate  
pH @ 20-25°C: 6.8 +/- 0.2

### **Microbial Load Evaluation:**

Manufactured lot has been statistically sampled. Samples were acceptable for microbial content and aesthetic value following a minimum of 7 days at 20-25C and 33-37C.

Signed



Quality Assurance Supervisor