

LOT

Certificate of Analysis

Product Name: Bruker Matrix HCCA, portioned

Part No.: 8255344

LOT No.: 7020120034

Use by¹⁾: 2022-03

Date of Production: 2020-04-14

Storage Conditions: Store at 2°C to 8°C on arrival.

Test:

Specification

Results

1. Amount of HCCA content.

Sampling in accordance to DIN ISO
U=2.8 mg // L=2.2 mg

2. MALDI Biotyper identification
Classification Result Report lists
Passed

of a 1:20 diluted BTS.

Sampling in accordance to DIN ISO

2859-1.

AQL-value=1.0
Inspection level S4

L= 2.000 log(score) of the best match

3. Blank value analysis of HCCAportioned by MALDI Biotyper.
Sampling in accordance to DIN ISO
2859-1.

Classification Result Report lists "no reliable identification".

Passed

Passed

Inspection level S4 U= 1.699 log(score)

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 22-04-2020

Fahrenheitstraße 4
D-28359 Bremen
Tel: +49 (0) 421 - 22 05 - 0
Dr.M. Kostrzewa, SAFETY OFFICER

Not for use in diagnostic procedures. For research use only.