

# SAFE® CAP 96 CHANNEL DD

## Automated 96-Channel Capper/Decapper

- Fully-automated 96 channel premium capper/decapper
- Process a full rack of tubes
- Caps can be hold or placed on a cap rack
- Can be operated as an integrated and standalone device
- Touchscreen panel provides easy navigation through the system
- Capping/Decapping time respectively approx. 50 seconds
- Latest cap driver technology: Cap drivers with spring-loaded metal balls optimally adapt to the caps and permit 100% leak-tight and material-saving sealing of the tubes without tilting
- Innovative and low-maintenance rack-to-head technology: It is not the complex driver head moves to the rack, but the rack is lifted to the head
- Automated height detection
- Automated rack barcode recognition
- Very small footprint 408 x 234 mm
- Height 388 mm

### Item Number

SAFE® 96 XT:

SAFE® 96 IT:

CDC-96CH-XT-DD

CDC-96CH-IT-DD

### Order Information



SAFE® CAP 96 CHANNEL DD