

# CERTIFICATE OF ANALYSIS

## POSITIVE CONTROL KIT

Rev. 1.0 - 2015.07.01

<b>PRODUCT CODE: SI 190.911</b>	<b>NAME: POSITIVE CONTROL KIT</b>
<b>BATCH N.: PCK16002</b>	
<b>MANUFACTURING DATE: 30/03/2016</b>	<b>EXPIRY DATE: 30/05/2017</b>

◆ *Reconstituting and Diluting Solutions*

Eugonic broth

Batch N.: FLA16001

Expire date:  
19/02/2018

	Physical parameters *	Decontamination at 37 °C after 72 h	Growth inter batch	Growth intra batch
<b>Expected</b>	Compliant	No growth	CV%≤35%	CV%≤35%
<b>Found</b>	OK	OK	7.4%	<5%

\* Volume, transparency, colour

◆ *Lyophilized Microorganism*

Batch N: 15003

Expire date: 01/12/2017

	Growth inter batch (CV%≤35%)	Physical check**
<b>Expected</b>	Compliant	Compliant
<b>Found</b>	OK	OK

\*\* Aspect

**STORAGE CONDITIONS: 2°C – 8 °C**

**Date of Approval: 04/04/2016**

Responsible for batch release


