

### Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,  
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella H g,m**

Catalogue No. (Package size) TR 1408 (1 ml)

Lot No. **54 05 13**

Expiry date **2015-06**

Intended use Enteroclon Anti-Salmonella H g,m records all Salmonella strains, as long as they contain in an antigen combination the H g and the H m antigens (e.g. f,g; f,g,s; f,g,t; g,s,t; g,q; g,t; g,z<sub>51</sub> or g,m; g,m,s,t; m,t) in slide agglutination.

Specificity No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2..8 °C

Berlin,

01.07.13

  
Head of Quality Control