

Certificate of Quality Control

| | |
|---------------------------------|---|
| Manufacturer | SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany |
| Product | Enteroclon Anti-Salmonella H g,m |
| Catalogue No. (Package size) | TR 1408 (1 ml) |
| Lot No. | 63 02 14 |
| Expiry date | 2015-09 |
| Intended use | Enteroclon Anti-Salmonella H g,m records all Salmonella strains, as long as they contain in an antigen combination the H g and the H m antigens (e.g. f,g; f,g,s; f,g,t; g,s,t; g,q; g,t; g,z ₅₁ or g,m; g,m,s,t; m,t) in slide agglutination. |
| Specificity | No agglutination with heterologous Salmonella H antigens and Salmonella O antigens. |
| Preservatives | Sodium azide 0.9 mg/ ml |
| Storage | 2...8 °C |

Berlin,

24.03.14



Head of Quality Control