

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:g,m</b>
Catalogue no.	TR 1408
Package size	1 ml
<b>Lot no.</b>	<b>8670917</b>
<b>Expiry date</b>	<b>2019-09</b>
Intended use	Anti-Salmonella H:g,m records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:g and/or H:m antigens (e. g. H:f,g; H:f,g,s; H:f,g,t; H:g,p; H:g,p,u; H:g,q; H:g,s,t; H:g,t; H:g,z <sub>51</sub> or H:g,m; H:g,m,s,t; H:m,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 26. 10. 20 17  
 place, date

J. Heu den  
 signature Head of Quality Control