

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:g,m
Catalogue no.	TR 1408
Package size	1 ml
Lot no.	7890919
Expiry date	2021-09
Intended use	Anti-Salmonella H:g,m records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:g and/or H:m antigens (e. g. H:f,g; H:f,g,s; H:f,g,t; H:g,p; H:g,p,u; H:g,q; H:g,s,t; H:g,t; H:g,z ₅₁ or H:g,m; H:g,m,s,t; H:m,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 22.10.19
place, date

7. Landerj
signature Head of Quality Control