

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:g,m
Catalogue no.	TR1408
Package size	1 ml
Lot no.	8221020
Expiry date	2022-09
Intended use	Anti-Salmonella H:g,m records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:g and/or H:m antigens (e. g. H:f,g; H:f,g,s; H:f,g,t; H:g,p; H:g,p,u; H:g,q; H:g,s,t; H:g,t; H:g,z ₅₁ or H:g,m; H:g,m,s,t; H:m,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 04.11.2020
place, date

J. Lendy
signature Head of Quality Control