

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485 sifin diagnostics gmbh

Berliner Allee 317-321 13088 Berlin, Germany

T: +49 30 927 030-0 F: +49 30 927 030-30 info@sifin.de www.sifin.de

**Product** 

Anti-Coli O 124:K 72

Catalogue no.

TR 2210

Package size

1 ml

Lot no.

22 01 15

**Expiry date** 

2020-01

Intended use

Anti-Coli O 124:K 72 records all E. coli strains that contain the antigen O 124:K 72 in slide agglutination. In the Widal reaction the serum reaches a titre of 1:1600 (tube and MTP) with the homologous E. coli strains; deviations

of ±1 titre step are acceptable.

Specificity

Since unspecific co-agglutinations with other E. coli strains often occur in the slide agglutination test, a positive result in the slide test must always be confirmed by examination of the boiled E. coli isolate in the Widal test.

Preservatives

Sodium azide 0.9 mg/ ml

Storage

2 °C...8 °C

Berlin, 04,06.15

place, date

signature Head of Quality Control