


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Serum Anti-Salmonella H m
Catalogue No. (Package size)	TS 1413 (1 ml)
Lot No.	126 12 10
Expiry date	2015-12
Intended use	Test Serum Anti-Salmonella H m records all Salmonella strains, as long as they contain in an antigen combination the H m antigen (e.g. H g,m; H g,m,s,t; H m,t) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

06.01.11


Head of Quality Control