

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test serum Anti-Yersinia enterocolitica O 3
Catalogue No. (Package size)	TS 1701 (1 ml)
Lot No.	82 05 12
Expiry date	2017-05
Intended use	Records all Yersinia enterocolitica strains in the slide agglutination as long as they contain the O 3 antigen
Specificity	No agglutination with heterologous Yersinia enterocolitica O antigens and selected representatives of other Enterobacteriaceae.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 03.07.12



Head of Quality Control