


Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Yersinia enterocolitica O 3
Catalogue no.	TS 1703
Package size	1 ml
Lot no.	176 08 15
Expiry date	2020-08
Intended use	Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 3 in slide agglutination.
Specificity	No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 19.10.15
place, date


signature Head of Quality Control