

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Yersinia enterocolitica O 8
Catalogue no.	TS 1705
Package size	1 ml
Lot no.	2600318
Expiry date	2023-03
Intended use	Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 8 in slide agglutination.
Specificity	No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 25.05.2018

 place, date

T. Schlüter

 signature Head of Quality Control