

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test serum Anti-Yersinia enterocolitica O 27**

Catalogue No. (Package size) TS 1706 (1 ml)

Lot No. **74 06 13**

Expiry date **2016-06**

Intended use Records all Yersinia enterocolitica strains in the slide agglutination as long as they contain the O 27 antigen in addition to the O 5 antigen.

Specificity No agglutination with heterologous Yersinia enterocolitica O antigens and selected representatives of other Enterobacteriaceae.

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin, 25.06.13



Head of Quality Control