

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Yersinia enterocolitica O 27</b>
Catalogue no.	TS 1706
Package size	1 ml
<b>Lot no.</b>	<b>66 03 15</b>
<b>Expiry date</b>	<b>2018-03</b>
Intended use	Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 27 in addition to antigen O 5 in slide agglutination.
Specificity	No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
place, date

17.04.15

signature Head of Quality Control

