

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Shigella flexneri type 6
Catalogue no.	TS 2008
Package size	1 ml
Lot no.	119 06 16
Expiry date	2021-06
Intended use	Anti-Shigella flexneri type 6 records all <i>Shigella flexneri</i> strains, as long as they contain the serotype 6 in slide agglutination.
Specificity	No agglutination with serotypes of <i>Shigella flexneri</i> types 1 to 5 and <i>Shigella sonnei</i> (phase I and phase II).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2 °C...8 °C

Berlin,
place, date

14.07.2016


signature Head of Quality Control