


Nutrient Broth N°2

N° Cat: 1251

N° Batch: 301314

Re-test Date: 1/2026



Formula in g/l

Casein peptone	10	Beef extract	10
Sodium chloride	5		

Typical formula g/L * Adjusted and/or supplemented as required to meet performance criteria.

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Batch pH
w/o rests	Fine powder	Beige	Amber, slightly opalescent	7,5 ± 0,2	7,3

Microbiological test

Incubation conditions: (35±2 °C / 18-48 h)

Microorganisms	Specification
Staphylococcus epidermidis ATCC 12228	Good growth
Klebsiella aerogenes ATCC 13048	Good growth
Streptococcus pyogenes ATCC 19615	Good growth
Escherichia coli ATCC 25922	Good growth
Salmonella typhi ATCC 6539	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.