

CERTIFICATE OF ANALYSIS

Product Name: SAFRANIN	Lot Number: K34024
Product Code: PL.7012	Expiration Date: 2023 07 31

Additional Information

Description: Deep red solution with slight odour.	Concentrate Lot number: K33957 Production date: 2022/07/26	Storage conditions: 15-25°C Precautions: refer to SDS
---------------------------------------------------	---------------------------------------------------------------	----------------------------------------------------------


Quality Control Testing

Challenged with test organisms:
Escherichia coli NCTC 12241
Oxford Staphylococcus aureus NCTC 12981

Test results:
 Conforms to performance specifications when used in standard Grams stain method, giving crisp staining of Gram positive and Gram negative bacteria with a clean counterstained background.

Comments:
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release

2022/07/29	
Date	Quality Control