

SAFE® READ Single Tube

2D Code Single Tube Reader

- Fast scanning of 2D tubes
- Compatible with SAFE® and other 2D tubes
- Keyboard Wedge setting allows data entry directly into your application
- Strong mineral glass window
- Decoding software with free lifetime upgrades
- Operating system on host PC = Windows 7, 8, 10
- Cryoprotection option
- Dimensions: 60mm x 90mm x 60mm
- 2 years warranty

Item Number

DMST

Order Information



SAFE® READ Single Tube