

Certificate of Analysis

Product Name:

Bruker Bacterial Test Standard

Part No.

8255343

LOT No.:

7030220016

Use by¹⁾:

2022-09-30

Date of Production:

2020-10-21

QC release date:

2020-10-26

Storage Conditions:

Store at -18 °C or below immediately on arrival.

<u>Test</u>	Specification	LOT Results
Verification of the calibration points by MALDI-TOF MS	Within a mass range of \pm 900ppm	Passed
Relative intensity of reference proteins No. 3 by MALDI-TOF MS	0.20 ± 0.05	Passed
Relative intensity of reference proteins No. 4 by MALDI-TOF MS	0.22 ± 0.05	Passed
MALDI Biotyper identification	Classification Result Report lists "Escherichia coli" as best match.	Passed
	log(score) Value ≥ 2.2 of the best match	Passed

Remarks:

1) BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker BTS can be stored for up to 5 months at -18°C or below. See instruction for use at www.bruker.com/care.

Bruker Daltonik GmbH

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

D-28359 Bremen Tel: +49(0) 421 + 22 05 0

Fahrenheitstraße

Dr. Markus Kostrzewa, SAFETY OFFICER