

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Yersinia enterocolitica O 3
Catalogue no.	TS1701
Package size	1 ml
Lot no.	7531022
Expiry date	2027-09
Intended use	Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 3 in slide agglutination.
Specificity	No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

**This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.**