
Biosan SIA, Ratsupites iela 7 k-2, Riga, Latvia, LV-1067 E-mail: marketing@biosan.lv, www.biosan.lv Phone:+37167426137



BN-13/75, adapter set

DESCRIPTION

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/75 for vacutainers 2-5 ml



BS-010208-PK

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-13/75 for vacutainers 2-5 ml





LMC-3000 Without rotor Laboratory Centrifuge

LMC-3000 is a modern lowspeed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.



LMC-12 Without rotor Laboratory Centrifuge

LMC-3000 is a modern lowspeed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

read more

Features: ...



LMC-4200R Without rotor Laboratory Refrigerated Centrifuge

Laboratory bench-top centrifuge with refrigeration LMC-4200R provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "coldshelf" is a gold standard for enzymologists and cell biologists because it ...





LMC-56 Without rotor Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

<u>read more</u>



R-12/10 BS-010208-BK Rotor

Swing-out rotor for round bottom centrifuge tube.

read more

Features: ...



Ri-12/10 BS-010221-MK Rotor

Swing-out rotor for round bottom centrifuge tube.