

## BN-16/100, adapter set

### DESCRIPTION

For rotors R-12/10, Ri-12/10 optional adapter set (12 pcs) - BN-16/100 for vacutainers 8-10 ml



### CAT. NUMBER

**BS-010208-RK**

For rotors R-12/10, Ri-12/10 optional adapter set  
(12 pcs) - BN-16/100 for vacutainers 8-10 ml



**LMC-3000**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-12**  
Without rotor  
Laboratory Centrifuge

**LMC-3000** is a modern low-speed bench-top centrifuge designed for operation with microtest plates and centrifuge tubes up to 50 ml. This device is widely used in biomedical profile laboratories.

Features: ...

[read more](#)



**LMC-4200R**  
Without rotor  
Laboratory Refrigerated Centrifuge

Laboratory bench-top centrifuge with refrigeration **LMC-4200R** provides temperature control of biomaterial during centrifugation. Temperature control of the so-called "cold-shelf" is a gold standard for enzymologists and cell biologists because it ...

[read more](#)



**LMC-56**  
Without rotor  
Laboratory Centrifuge

**LMC-56** is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



**R-12/10**  
BS-010208-BK  
Rotor

Swing-out rotor for round bottom centrifuge tube.



**Ri-12/10**  
BS-010221-MK  
Rotor

Swing-out rotor for round bottom centrifuge tube.