

Microblot-Array

Multiplex diagnostics
in microtiter plate format

OVERVIEW




Multiplex testing



ELISA-compatible



Efficiency



Quantitative results



Easy-to-use software




Find out more

Multiplex diagnostics

Microblot-Array is an immunoblot array in the well of a microtiter plate format designed for efficient multiplex diagnostics. It brings ELISA's high throughput to BLOT testing. The comprehensive evaluation is ensured

by using the TestLine Analytics Software in combination with the Microblot-Array Reader, enabling complex image analysis, including results evaluation and connectivity to LIS.

So much to discover – in a single well.



Innovation is part of our portfolio.

Panel	Parameters included	Availability
Atypical respirations	Mycoplasma pneumoniae, Bordetella, Chlamydia pneumoniae, Legionella pneumophila	In development
Gastrointestinal diseases	Helicobacter pylori, Yersinia enterocolitica, Autoimmune gastrointestinal diseases	In development
Autoimmune neurological diseases	Paraneoplastic syndrome, Limbic encephalitis	In development
Herpetic infections	EBV, CMV, HSV, VZV, HHV-6	In development
ANCA	MPO, PR3, GBM, Lactoferin, Lysozyme, BPI, Cathepsin, Elastase	In development
Vector transmitted infections	Rickettsia, Babesia, Anaplasma, Neoehrlichia, TBEV, Borrelia, Francisella, Coxiella burnetti	In development
Tropical diseases	Dengue, Chikungunya, Zika, West Nile Fever, Plasmodium, Rickettsia, Leptospira	In development

The composition of the panels can be changed

4 steps, countless options.

1

Unbox


Ready to use kit,
no dilution needed



2

Compose

Up to 96 various tests in one plate
thanks to breakable wells



3

Process

Automated or manual processing
with standard ELISA routine



4

Evaluate

Comprehensive evaluation with
dedicated MBA Reader and Software





FOLLOW US

BIOVENDOR.GROUP



Diagnostic kits are intended
for professional use in the laboratory.

Portfolio. Unique composition, comprehensive diagnostics.

AUTOIMMUNITY

Code	Products	No. of tests
ANApMA96	Microblot-Array ANA plus*	96
LKMMA48	Microblot-Array Liver profile	48
AIGAMA48	Microblot-Array Autoimmune gastroenteritis panel IgA	48
AIGGMA48	Microblot-Array Autoimmune gastroenteritis panel IgG	48

*Check availability in your country.

INFECTIOUS SEROLOGY

Code	Products	No. of tests
BpAMA48	Microblot-Array Bordetella IgA	48
BpGMA48	Microblot-Array Bordetella IgG	48
BpMMA48	Microblot-Array Bordetella IgM	48
BGMA096	Microblot-Array Borrelia IgG	96
BMMA096	Microblot-Array Borrelia IgM	96
CoVAMA96	Microblot-Array COVID-19 IgA	96
CoVGMA96	Microblot-Array COVID-19 IgG	96
CoVMMA96	Microblot-Array COVID-19 IgM	96
CMGMA48	Microblot-Array CMV IgG	48
CMMMA48	Microblot-Array CMV IgM	48
EBAMA96	Microblot-Array EBV IgA	96
EBGMA96	Microblot-Array EBV IgG	96
EBMMA96	Microblot-Array EBV IgM	96
HpAMA48	Microblot-Array Helicobacter IgA	48
HpGMA48	Microblot-Array Helicobacter IgG	48
HSGMA48	Microblot-Array HSV 1+2 IgG	48
HSMMA48	Microblot-Array HSV 1+2 IgM	48
CAMA096	Microblot-Array Chlamydia IgA	96
CGMA096	Microblot-Array Chlamydia IgG	96
MyAMA48	Microblot-Array Mycoplasma IgA	48
MyGMA48	Microblot-Array Mycoplasma IgG	48
MyMMA48	Microblot-Array Mycoplasma IgM	48
YAMA048	Microblot-Array Yersinia IgA	48
YGMA048	Microblot-Array Yersinia IgG	48

PRODUCER

TestLine Clinical Diagnostics s.r.o.

Křižíkova 188/68

612 00 Brno

Czech Republic