

CERTIFICATE OF ANALYSIS

Product Name:	Staph Latex Kit		Lot Number:			C22210	
Product Code:	PL.080B			Expir	ation Date:	2025-04-30	
Product Code	Quantity	Compone	nts				
PL.083B	2 x 2.5 ml	Staph Late	Staph Latex Reagent				
PL.085B	1 x 2.5 ml	Negative C	Negative Control Latex Reagent				
	200	Mixing Stic	ks				
PL.092-20	20 cards	Test Cards	Cards				
Additional Inform							
Storage Conditions: 2-8°C							
Instructions for Use: Refer to pro-lab.com/resources for current version							
Quality Control Testing							
Criteria						Result	
Strong agglutination of Staph Latex Reagent when tested against Staphylococcus aureus.						PASS	
No agglutination of Negative Control when tested against test isolates.						PASS	
Release							
2023-11-22			L	ava	lu		
Date				Quality Control			