



www.pro-lab.co.uk

Certificate of analysis

Product Name: MALACHITE GREEN	Lot number: K38900
Product Code: PL.8007	Expiry date: 2025 09 30

Additional Information		
Description: Dark green solution.	Concentrate Lot number: K38792 Production date: 2024/09/09 Instructions for Use:	Storage conditions: 15-25°C Precautions: refer to SDS Go to pro-lab.co.uk for current version.

Quality Control Testing
<p>Challenged with test organisms: <i>Mycobacterium scrofulaceum</i> NCTC 10803 <i>Escherichia coli</i> NCTC 12241</p> <p>Test results: Conforms to performance specifications when used in standard Ziehl Neelson stain method, giving crisp red staining of Acid/AlcoholFast bacilli with a clean counterstained background.</p> <p>Comments: We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.</p>

Release	
2024-09-10	R. Grover
Date	Signed

T: 0151 353 1613
3 Bassendale Road, Bromborough,
Wirral, CH62 3QL, UK.
uksupport@pro-lab.com