



www.pro-lab.co.uk

Certificate of analysis

Product Name: LUGOLS IODINE	Lot number: K38710
Product Code: PL.8010	Expiry date: 2025 08 31

Additional Information		
Description: Dark brown solution.	Concentrate Lot number: K38670 Production date: 2024/08/08 Instructions for Use:	Storage conditions: 15-25°C Precautions: refer to SDS Go to pro-lab.co.uk for current version.

Quality Control Testing
<p>Challenged with test organisms: <i>Escherichia coli</i> NCTC 12241 Oxford <i>Staphylococcus aureus</i> NCTC 12981</p> <p>Test results: Conforms to performance specifications when used in standard Grams stain method, giving crisp staining of Gram positive and Gram negative bacteria with a clean counterstained background.</p> <p>Comments: We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.</p>

Release	
2024-08-13	R. Grover
Date	Signed

T: 0151 353 1613
3 Bassendale Road, Bromborough,
Wirral, CH62 3QL, UK.
uksupport@pro-lab.com