

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Yersinia enterocolitica O 8
Catalogue no.	TS1705
Package size	1 ml
Lot no.	3740524
Expiry date	2029-04
Intended use	Records all <i>Yersinia enterocolitica</i> strains as long as they contain antigen O 8 in slide agglutination.
Specificity	No agglutination with heterologous <i>Yersinia enterocolitica</i> O antigens and selected representatives of other <i>Enterobacteriaceae</i> .
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

This document has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.