

Certificate of Analysis

PRODUCT DESCRIPTION: NADAL® hCG Pregnancy Test

PRODUCT	LOT	LOT AMOUNT	EXPIRATION	PHYSICAL
CODE	NUMBER	LOT AMOUNT	DATE	APPEARANCE
152002	175845	24000PCS	2024-01	Conforms

PRODUCT TEST RESULTS

TEST		STANDARDS		QC CRITERION	RESULTS	PASS (Y/N)
Specificity:	300 mIU/ml	LH (2nd IS)	in Buffer	Negative	Negative	YES
Sensitivity 510(k) claimed:	25mIU/ml	HCG(4th IPR)	in Buffer	Positive	Positive	YES
Prozone:	600 IU/ml	HCG(4th IPR)	in Buffer	Positive	Positive	YES
For information only:	25mIU/mI	HCG(4th IPR)	in Urine	Positive	Positive	YES

PRODUCT MANUFACTURING:

No radioactive material of any kind is utilized in the product, or in the manufacturing of the product. The product specified above does not contain any virus, reagent by-product of the same, or metabolic by-product of Hepatitis A, B, C, or D.

If this product is shipped as bulk, nal von minden GmbH assumes no responsibility for the packaging and label claims.

Accepted:

Sandra von Minden

Qualitätsmanagement nal von minden GmbH nal von minden Gmb Carl-Zeiss-Straße 12 47445 Moers Tel: +49 2841 99820-0 Fax: +49 2841 99820-1

von minden info@nal-vonminden.de

Moers, 2021-08-03

nal von minden GmbH

Carl-Zeiss-Straße 12 47445 Moers Deutschland

Friedenstraße 32 93053 Regensburg Deutschland

Tel.: +49 941 29010-0 Fax: +49 941 29010-50

Geschäftsführer: Sandra von Minden info@nal-vonminden.com Roland Meißner www.nal-vonminden.com Thomas Zander

Registergericht Kleve: HRB 5679

Sparkasse Regensburg Deutsche Bank Steuer-Nr. 244/133/00130 Kto: 840 017 073 UST-ID-Nr. DE 189 016 086 BLZ: 750 500 00

Kto: 056 998 800 BLZ: 350 700 24

