

Certificate of Analysis

PRODUCT DESCRIPTION: NADAL® COVID-19 Ag plus Test

PRODUCT	LOT	EXPIRATION	PHYSICAL
CODE	NUMBER	DATE	APPEARANCE
243104N-20	SE2011197	2022-12-04	Conforms

PRODUCT TEST RESULTS

Inspection item	Acceptance Criteria		Test Result	
Appearance	The test board has a flat appearance, is firmly pasted, has no smudges, and is closely attached. Extraction solution appearance: No turbidity, no precipitation, no leakage, no stains on the appearance. The kit is not damaged or crushed, and all information is clearly marked.		PASS	
Inspection item	Reference samples	The test results compare with the color grade		Test result
		Evaluation criteria of C line	Evaluation criteria of T line	
Positive reference samples	P1	≥G9	G3-G4	PASS
	P2		G4-G5	PASS
	P3		G5-G6	PASS
	P4		G6-G7	PASS
	P5		G7-G8	PASS
	P6		G8-G9	PASS
	P7		G9-G10	PASS
	P8		>G10	PASS
Negative reference samples	N1-N20	≥G9	G1-G2	PASS
Minimum detection limit reference samples	S1	≥G9	G6-G7	PASS
	S2		G5-G6	PASS
	S3		G4-G5	PASS
	S4		G3-G4	PASS
	S5		G3	PASS
	S6		G1-G2	PASS
Repeatability reference samples	R1	≥G9	G3-G4	PASS
	R2		G5-G6	PASS
Conclusion	PASS			

nal von minden GmbH

Carl-Zeiss-Straße 12
47445 Moers
Deutschland

Friedenstraße 32
93053 Regensburg
Deutschland

Tel.: +49 941 29010-0
Fax: +49 941 29010-50

info@nal-vonminden.com
www.nal-vonminden.com

Geschäftsführer:
Sandra von Minden
Roland Meißner
Thomas Zander

Registergericht Kleve:
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

Sparkasse Regensburg
Kto: 840 017 073
BLZ: 750 500 00

Deutsche Bank
Kto: 056 998 800
BLZ: 350 700 24

Schnelltests

Labordiagnostika

Laborservice

Beratung & Service

PRODUCT MANUFACTURING:

No radioactive material of any kind is utilized in the product, or in the manufacturing of the product. The product specified above does not contain any virus, reagent by-product of the same, or metabolic by-product of Hepatitis A, B, C, or D.

If this product is shipped as bulk, nal von minden GmbH assumes no responsibility for the packaging and label claims.

Accepted:

Sandra von Minden

Qualitätsmanagement
nal von minden GmbH



nal von minden GmbH
Carl-Zeiss-Straße 12
47445 Moers
Tel: +49 2841 99820-0
Fax: +49 2841 99820-1
info@nal-vonminden.de



nal
von minden

Moers, 2021-06-25

nal von minden GmbH

Carl-Zeiss-Straße 12
47445 Moers
Deutschland

Friedenstraße 32
93053 Regensburg
Deutschland

Tel.: +49 941 29010-0
Fax: +49 941 29010-50

info@nal-vonminden.com
www.nal-vonminden.com

Geschäftsführer:
Sandra von Minden
Roland Meißner
Thomas Zander

Registergericht Kleve:
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

Sparkasse Regensburg
Kto: 840 017 073
BLZ: 750 500 00

Deutsche Bank
Kto: 056 998 800
BLZ: 350 700 24

Schnelltests

Labordiagnostika

Laborservice

Beratung & Service