

Certificate of Quality

Product Name: VERSATREK REDOX 2, EZ DRAW 40ML

Lot Number: 493718

Product Number: 7107-44 Expiration Date: 2015-04-28

(YYYY-MM-DD)

This product has been manufactured, processed and packaged in accordance with Quality Systems Regulation, 21 CFR Part 820. Representative samples were tested per Remel Quality Control specifications and were found to meet performance criteria for this product.

Organism Challenge

ATCC

Performance Criteria

Bacteroides fragilis

25285

Positive within 96 Hours

Streptococcus pneumoniae

6305

Positive within 96 Hours

Physical Characteristics:

Appearance:

Light Amber, Clear Broth w/ No Pink-Red Discoloration or Significant Precipitate

pH @ 20-25°C: 6.9 +/- 0.2

Microbial Load Evaluation:

Manufactured lot has been statistically sampled. Samples were acceptable for microbial content and aesthetic value following a minimum of 72 hours at 33-37°C.

Signed

Susie Waddle

Quality Assurance Supervisor