

Certificate of Quality

Product Name: VERSATREK MYCO PVNA, 5/BX
Lot Number: 982177

Product Number: 7113-42
Expiration Date: 2022-02-28
(YYYY-MM-DD)

This product has been manufactured, processed and packaged in accordance with Quality Systems Regulation, 21 CFR Part 820. Representative samples were tested per Remel Inc., a part of Thermo Fisher Scientific Quality Control specifications and were found to meet performance criteria for this product.

Performance:

Solubility:	Soluble When Reconstituted	
Organism:	MIC Range:	Antibiotic:
Pseudomonas aeruginosa ATCC 27853	1.56 - 5 Units/mL	Polymyxin B
Staphylococcus aureus CDC 2081	0.75 - 1.5 µg/mL	Vancomycin
Proteus vulgaris ATCC 13315	2.5 - 5.0 µg/mL	Nalidixic Acid
Candida albicans ATCC 10231	0.1 - 0.2 µg/mL	Amphotericin B

Physical Characteristics:

Appearance: Yellow, Clear Liquid

Microbial Load Evaluation:

Samples were acceptable for microbial content following a minimum of 14 days at 33-37C.

Signed



Quality Assurance Supervisor