QUALITY CONTROL CERTIFICATE

AMPHOTERICIN B 10 µg Neo-Sensitabs

REF: 81012

LOT: 1901-2

Date of analysis: 2019-05

Date of expiry: 2022-04

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 90028 Candida parapsolis ATCC 22019 Candida krusei ATCC 6258

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Microbiological quality: Complies.

Manager Quality Assurance

Rosco Diagnostica A/S

Peter Samuelsen

Manager Quality Assurance